

Automated Microfluidic Platform

1.0

Generated by Doxygen 1.8.15

1 Namespace Index	1
1.1 Namespace List	1
2 Hierarchical Index	3
2.1 Class Hierarchy	3
3 Class Index	5
3.1 Class List	5
4 File Index	7
4.1 File List	7
5 Namespace Documentation	9
5.1 Collection Namespace Reference	9
5.1.1 Detailed Description	9
5.2 fakeSerial Namespace Reference	9
5.2.1 Detailed Description	9
5.3 gui Namespace Reference	9
5.3.1 Detailed Description	10
5.3.2 Function Documentation	10
5.3.2.1 check_input()	10
5.3.2.2 combine_funcs()	10
5.3.2.3 message_prompt()	10
5.3.3 Variable Documentation	11
5.3.3.1 app	11
5.3.3.2 LARGE_FONT	11
5.4 Mswitch Namespace Reference	11
5.4.1 Detailed Description	11
5.5 Pump Namespace Reference	11
5.5.1 Detailed Description	11
5.5.2 Function Documentation	12
5.5.2.1 get_ports()	12
5.6 TwoSwitch Namespace Reference	12
5.6.1 Detailed Description	12
6 Class Documentation	13
6.1 gui.App Class Reference	13
6.1.1 Detailed Description	14
6.1.2 Constructor & Destructor Documentation	14
6.1.2.1 __init__()	14
6.1.3 Member Function Documentation	14
6.1.3.1 done_and_load()	14
6.1.3.2 get_comport()	15
6.1.3.3 get_ports()	15

6.1.3.4 message_window()	15
6.1.3.5 negative_run_status()	15
6.1.3.6 positive_run_status()	16
6.1.3.7 set_defaults()	16
6.1.3.8 show_frame()	16
6.1.4 Member Data Documentation	16
6.1.4.1 appHighlightFont	16
6.1.4.2 buttonFont	16
6.1.4.3 default_diameter_index	17
6.1.4.4 frames	17
6.1.4.5 hasColl	17
6.1.4.6 hasMani	17
6.1.4.7 hasPump	17
6.1.4.8 max_flowrate_1	17
6.1.4.9 max_flowrate_2	17
6.1.4.10 max_flowrate_3	18
6.1.4.11 max_flowrate_list	18
6.1.4.12 my2Switch	18
6.1.4.13 myColl	18
6.1.4.14 myMani	18
6.1.4.15 myPump	18
6.1.4.16 pageFont	18
6.1.4.17 port2Switch	18
6.1.4.18 portColl	19
6.1.4.19 portMani	19
6.1.4.20 portPump	19
6.1.4.21 ports	19
6.1.4.22 should_be_running	19
6.1.4.23 well_working_volume	19
6.2 guiAutomaticPage Class Reference	19
6.2.1 Constructor & Destructor Documentation	20
6.2.1.1 __init__()	20
6.2.2 Member Function Documentation	20
6.2.2.1 add_step()	21
6.2.2.2 check_channel()	21
6.2.2.3 clear_recipe()	21
6.2.2.4 create_step()	21
6.2.2.5 delete_step()	21
6.2.2.6 file_load()	21
6.2.2.7 file_save()	22
6.2.2.8 run_recipe()	22
6.2.3 Member Data Documentation	22

6.2.3.1 channels	22
6.2.3.2 clearbutton	22
6.2.3.3 deletestep	22
6.2.3.4 deletestepcombo	22
6.2.3.5 ejectbutton	23
6.2.3.6 headinglabel	23
6.2.3.7 loadbutton	23
6.2.3.8 mainlabel	23
6.2.3.9 manualbutton	23
6.2.3.10 newstep	23
6.2.3.11 number_deleted	23
6.2.3.12 numberlist	23
6.2.3.13 reservoirs	24
6.2.3.14 runbutton	24
6.2.3.15 savebutton	24
6.2.3.16 settingsbutton	24
6.2.3.17 steplist	24
6.2.3.18 stepnumbers	24
6.2.3.19 todelete	24
6.2.3.20 twoSwitchValues	25
6.3 Collection.Collectador Class Reference	25
6.3.1 Detailed Description	25
6.3.2 Constructor & Destructor Documentation	26
6.3.2.1 __init__()	26
6.3.3 Member Function Documentation	26
6.3.3.1 chop_return()	26
6.3.3.2 eject()	26
6.3.3.3 get_info()	27
6.3.3.4 last_site()	27
6.3.3.5 moveOneDown()	27
6.3.3.6 moveOneLeft()	27
6.3.3.7 moveOneRight()	27
6.3.3.8 moveOneUp()	28
6.3.3.9 next_site()	28
6.3.3.10 reset()	28
6.3.3.11 send()	28
6.3.3.12 serial_connect()	29
6.3.3.13 setOrigin()	29
6.3.3.14 toggle_pattern()	29
6.3.4 Member Data Documentation	29
6.3.4.1 baud	29
6.3.4.2 connected	29

6.3.4.3 currentPattern	30
6.3.4.4 port	30
6.3.4.5 position	30
6.3.4.6 ser	30
6.3.4.7 serialConnected	30
6.3.4.8 uniqueID	30
6.3.4.9 usingSnakePattern	30
6.4 gui.ManualPage Class Reference	31
6.4.1 Constructor & Destructor Documentation	32
6.4.1.1 __init__()	32
6.4.2 Member Function Documentation	32
6.4.2.1 check_all_manual()	32
6.4.2.2 check_all_manual_and_sampling()	32
6.4.2.3 check_channel()	33
6.4.2.4 create_channels()	33
6.4.2.5 done()	33
6.4.2.6 eject()	33
6.4.2.7 moveOneDown()	33
6.4.2.8 moveOneLeft()	34
6.4.2.9 moveOneRight()	34
6.4.2.10 moveOneUp()	34
6.4.2.11 next()	34
6.4.2.12 prev()	34
6.4.2.13 reset()	34
6.4.2.14 set_and_start()	35
6.4.2.15 set_and_start_all()	35
6.4.2.16 setOrigin()	35
6.4.2.17 start_and_collect()	35
6.4.2.18 stop()	35
6.4.2.19 stop_all()	36
6.4.2.20 toggle()	36
6.4.3 Member Data Documentation	36
6.4.3.1 autocollectlabel	36
6.4.3.2 autopagebutton	36
6.4.3.3 channels	36
6.4.3.4 channeltextlist	36
6.4.3.5 CollectionControlLabel	37
6.4.3.6 ejectbutton	37
6.4.3.7 flowrateentrylist	37
6.4.3.8 greaterlabelframe	37
6.4.3.9 labelframelist	37
6.4.3.10 lastwellbutton	37

6.4.3.11 mainlabel	37
6.4.3.12 nextwellbutton	37
6.4.3.13 patternLabel	38
6.4.3.14 positionLabel	38
6.4.3.15 reservoircombobox	38
6.4.3.16 reservoirlabel	38
6.4.3.17 reservoirs	38
6.4.3.18 reservoirvalues	38
6.4.3.19 resetbutton	38
6.4.3.20 sampleslabel	38
6.4.3.21 samplingrate	39
6.4.3.22 samplingrateentry	39
6.4.3.23 setOriginButton	39
6.4.3.24 settingsbutton	39
6.4.3.25 sites	39
6.4.3.26 sitesentry	39
6.4.3.27 startallbutton	39
6.4.3.28 startandcollectbutton	39
6.4.3.29 startbuttonlist	40
6.4.3.30 stopallbutton	40
6.4.3.31 stopbuttonlist	40
6.4.3.32 togglebutton	40
6.4.3.33 toplabels	40
6.4.3.34 twoSwitchComboBoxList	40
6.4.3.35 twoSwitchValues	40
6.5 Mswitch.MSwitch Class Reference	41
6.5.1 Detailed Description	41
6.5.2 Constructor & Destructor Documentation	41
6.5.2.1 __init__()	41
6.5.3 Member Function Documentation	41
6.5.3.1 chop_return()	42
6.5.3.2 get_info()	42
6.5.3.3 send()	42
6.5.3.4 serial_connect()	42
6.5.3.5 set_reservoir()	43
6.5.4 Member Data Documentation	43
6.5.4.1 baud	43
6.5.4.2 connected	43
6.5.4.3 port	43
6.5.4.4 res	43
6.5.4.5 ser	43
6.5.4.6 serialConnected	44

6.5.4.7 uniqueID	44
6.6 fakeSerial.Serial Class Reference	44
6.6.1 Detailed Description	45
6.6.2 Constructor & Destructor Documentation	45
6.6.2.1 __init__()	45
6.6.3 Member Function Documentation	45
6.6.3.1 __str__()	45
6.6.3.2 close()	46
6.6.3.3 flush()	46
6.6.3.4 flushInput()	46
6.6.3.5 flushOutput()	46
6.6.3.6 isOpen()	46
6.6.3.7 open()	47
6.6.3.8 read()	47
6.6.3.9 readline()	47
6.6.3.10 write()	47
6.6.4 Member Data Documentation	48
6.6.4.1 baudrate	48
6.6.4.2 bytesize	48
6.6.4.3 name	48
6.6.4.4 parity	48
6.6.4.5 port	48
6.6.4.6 rtcts	48
6.6.4.7 stopbits	48
6.6.4.8 timeout	49
6.6.4.9 xonxoff	49
6.7 gui.SettingsPage Class Reference	49
6.7.1 Detailed Description	50
6.7.2 Constructor & Destructor Documentation	50
6.7.2.1 __init__()	50
6.7.3 Member Function Documentation	50
6.7.3.1 calibrate()	50
6.7.3.2 update()	50
6.7.4 Member Data Documentation	50
6.7.4.1 autopagebutton	50
6.7.4.2 calibratebutton	51
6.7.4.3 calibratecombobox	51
6.7.4.4 calibrationlabel	51
6.7.4.5 channelcaliblabel	51
6.7.4.6 channels	51
6.7.4.7 channelvalues	51
6.7.4.8 diameterlabel	51

6.7.4.9 diameterlist	51
6.7.4.10 diametervalues	52
6.7.4.11 directionheading	52
6.7.4.12 directionlabellist	52
6.7.4.13 directions	52
6.7.4.14 directionslist	52
6.7.4.15 entercalibratedbutton	52
6.7.4.16 mainlabel	52
6.7.4.17 manualbutton	52
6.7.4.18 measuredVolume	53
6.7.4.19 measuredVolumeEntry	53
6.7.4.20 reservoirs	53
6.7.4.21 updatebutton	53
6.8 Pump.ThreePump Class Reference	53
6.8.1 Detailed Description	54
6.8.2 Constructor & Destructor Documentation	54
6.8.2.1 __init__()	54
6.8.3 Member Function Documentation	54
6.8.3.1 abort_calibration()	54
6.8.3.2 calibrate()	55
6.8.3.3 chop_return()	55
6.8.3.4 FormatVolume()	55
6.8.3.5 send()	55
6.8.3.6 send_return()	56
6.8.3.7 serial_connect()	56
6.8.3.8 setDefaults()	56
6.8.3.9 setDir()	56
6.8.3.10 setFlow()	57
6.8.3.11 setMeasuredVolume()	57
6.8.3.12 setTargetCalibrationVolume()	57
6.8.3.13 start()	57
6.8.3.14 start_all()	58
6.8.3.15 stop()	58
6.8.3.16 stop_all()	58
6.8.4 Member Data Documentation	58
6.8.4.1 baud	58
6.8.4.2 calibrationunit	58
6.8.4.3 calibrationvolume	59
6.8.4.4 connected	59
6.8.4.5 isOn	59
6.8.4.6 port	59
6.8.4.7 ser	59

6.8.4.8 serialConnected	59
6.9 TwoSwitch.TwoSwitch Class Reference	59
6.9.1 Detailed Description	60
6.9.2 Constructor & Destructor Documentation	60
6.9.2.1 __init__()	60
6.9.3 Member Function Documentation	60
6.9.3.1 chop_return()	60
6.9.3.2 getIdentifier()	61
6.9.3.3 send()	61
6.9.3.4 serial_connect()	61
6.9.3.5 setCollect()	61
6.9.3.6 setRecirculate()	62
6.9.4 Member Data Documentation	62
6.9.4.1 baud	62
6.9.4.2 connected	62
6.9.4.3 port	62
6.9.4.4 ser	62
6.9.4.5 serialConnected	62
6.9.4.6 uniqueID	63
6.9.4.7 willRecirculate	63
6.10 gui.WelcomePage Class Reference	63
6.10.1 Constructor & Destructor Documentation	63
6.10.1.1 __init__()	63
6.10.2 Member Function Documentation	64
6.10.2.1 flip()	64
6.10.3 Member Data Documentation	64
6.10.3.1 coll_label	64
6.10.3.2 collectorframe	64
6.10.3.3 donebutton	64
6.10.3.4 mainlabel	64
6.10.3.5 manilabel	64
6.10.3.6 perfusorframe	65
6.10.3.7 pumplabel	65
6.10.3.8 sublabelframe	65
7 File Documentation	67
7.1 ArduinoSnakePattern.cpp File Reference	67
7.1.1 Detailed Description	69
7.1.2 Function Documentation	69
7.1.2.1 backwardstep1()	69
7.1.2.2 backwardstep2()	69
7.1.2.3 eject()	69

7.1.2.4 forwardstep1()	69
7.1.2.5 forwardstep2()	70
7.1.2.6 loop()	70
7.1.2.7 moveDown()	70
7.1.2.8 moveGen()	70
7.1.2.9 moveLast()	70
7.1.2.10 moveLeft()	71
7.1.2.11 moveNext()	71
7.1.2.12 moveRight()	71
7.1.2.13 moveUp()	71
7.1.2.14 populateLocations()	71
7.1.2.15 reset()	72
7.1.2.16 set()	72
7.1.2.17 setup()	72
7.1.2.18 stepperB()	72
7.1.2.19 stepperT()	72
7.1.3 Variable Documentation	72
7.1.3.1 acceleration	72
7.1.3.2 AFMS	73
7.1.3.3 columns	73
7.1.3.4 delayRPM	73
7.1.3.5 limitSwitchB	73
7.1.3.6 limitSwitchT	73
7.1.3.7 locations	73
7.1.3.8 maximumSpeed	74
7.1.3.9 myMotorB	74
7.1.3.10 myMotorT	74
7.1.3.11 myPosition	74
7.1.3.12 numLocations	74
7.1.3.13 positions	74
7.1.3.14 rows	75
7.1.3.15 steppers	75
7.1.3.16 stepsPerRevolution	75
7.1.3.17 stepsPerWell	75
7.1.3.18 stepsToEject	75
7.1.3.19 TwoSwitchA	75
7.1.3.20 TwoSwitchB	76
7.1.3.21 TwoSwitchC	76
7.2 Collection.py File Reference	76
7.3 fakeSerial.py File Reference	76
7.4 gui.py File Reference	76
7.5 locationPopulate.cpp File Reference	77

7.5.1 Function Documentation	77
7.5.1.1 main()	77
7.5.1.2 populateLocations()	78
7.5.2 Variable Documentation	78
7.5.2.1 columns	78
7.5.2.2 numLocations	78
7.5.2.3 rows	78
7.5.2.4 stepsPerWell	78
7.5.2.5 stepsToEject	78
7.6 Mswitch.py File Reference	79
7.7 Pump.py File Reference	79
7.8 TwoSwitch.py File Reference	79
7.9 XYStage_LimitSwitch_TwoSwitch_Control.cpp File Reference	79
7.9.1 Function Documentation	81
7.9.1.1 backwardstep1()	81
7.9.1.2 backwardstep2()	81
7.9.1.3 eject()	82
7.9.1.4 forwardstep1()	82
7.9.1.5 forwardstep2()	82
7.9.1.6 loop()	82
7.9.1.7 moveDown()	82
7.9.1.8 moveGen()	82
7.9.1.9 moveLast()	83
7.9.1.10 moveLeft()	83
7.9.1.11 moveNext()	83
7.9.1.12 moveRight()	83
7.9.1.13 moveUp()	83
7.9.1.14 populateLocations()	84
7.9.1.15 reset()	84
7.9.1.16 set()	84
7.9.1.17 setup()	84
7.9.1.18 stepperB()	84
7.9.1.19 stepperT()	84
7.9.2 Variable Documentation	84
7.9.2.1 acceleration	85
7.9.2.2 AFMS	85
7.9.2.3 collectionPattern	85
7.9.2.4 columns	85
7.9.2.5 delayRPM	85
7.9.2.6 limitSwitchB	85
7.9.2.7 limitSwitchT	86
7.9.2.8 locations	86

7.9.2.9 maximumSpeed	86
7.9.2.10 myMotorB	86
7.9.2.11 myMotorT	86
7.9.2.12 myPosition	86
7.9.2.13 numLocations	87
7.9.2.14 positions	87
7.9.2.15 rows	87
7.9.2.16 steppers	87
7.9.2.17 stepsPerRevolution	87
7.9.2.18 stepsPerWell	87
7.9.2.19 stepsToEject	88
7.9.2.20 TwoSwitchA	88
7.9.2.21 TwoSwitchB	88
7.9.2.22 TwoSwitchC	88
Index	89

Chapter 1

Namespace Index

1.1 Namespace List

Here is a list of all namespaces with brief descriptions:

Collection	9
fakeSerial	9
gui	9
Mswitch	11
Pump	11
TwoSwitch	12

Chapter 2

Hierarchical Index

2.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

Collection.Collectador	25
Frame	
guiAutomaticPage	19
guiManualPage	31
guiSettingsPage	49
guiWelcomePage	63
Mswitch.MSwitch	41
fakeSerial.Serial	44
Pump.ThreePump	53
Tk	
gui.App	13
TwoSwitch.TwoSwitch	59

Chapter 3

Class Index

3.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

gui.App	13
guiAutomaticPage	19
Collection.Collectador	25
gui.ManualPage	31
Mswitch.MSwitch	41
fakeSerial.Serial	44
gui.SettingsPage	49
Pump.ThreePump	53
TwoSwitch.TwoSwitch	59
gui.WelcomePage	63

Chapter 4

File Index

4.1 File List

Here is a list of all files with brief descriptions:

ArduinoSnakePattern.cpp	
Arduino code for XYStage Movement and TwoSwitch control	67
Collection.py	76
fakeSerial.py	76
gui.py	76
locationPopulate.cpp	77
Mswitch.py	79
Pump.py	79
TwoSwitch.py	79
XYStage_LimitSwitch_TwoSwitch_Control.cpp	79

Chapter 5

Namespace Documentation

5.1 Collection Namespace Reference

Classes

- class [Collectador](#)

5.1.1 Detailed Description

Fluidic Handling Software
Ashutosh Agarwal Lab
University of Miami

by:
Liev Birman
Adiel Hernandez

5.2 fakeSerial Namespace Reference

Classes

- class [Serial](#)

5.2.1 Detailed Description

Written by:
D. Thiebaut

5.3 gui Namespace Reference

Classes

- class [App](#)
- class [AutomaticPage](#)
- class [ManualPage](#)
- class [SettingsPage](#)
- class [WelcomePage](#)

Functions

- def `combine_funcs` (*funcs)
- def `message_prompt` (message)
- def `check_input` (variable, desired_datatype, lower_limit=0, upper_limit=0)

Variables

- tuple `LARGE_FONT` = ("Verdana", 10)
- `app` = `App()`

5.3.1 Detailed Description

Fluidic Handling Software
Ashutosh Agarwal Lab
University of Miami

by:
Liev Birman
Adiel Hernandez

5.3.2 Function Documentation

5.3.2.1 `check_input()`

```
def gui.check_input (
    variable,
    desired_datatype,
    lower_limit = 0,
    upper_limit = 0 )
```

5.3.2.2 `combine_funcs()`

```
def gui.combine_funcs (
    * funcs )
```

Use this function for executing multiple functions simultaneously with Tkinter buttons

5.3.2.3 `message_prompt()`

```
def gui.message_prompt (
    message )
```

Creates a pop-up message of your choice

5.3.3 Variable Documentation

5.3.3.1 app

```
gui.app = App()
```

5.3.3.2 LARGE_FONT

```
tuple gui.LARGE_FONT = ("Verdana", 10)
```

5.4 Mswitch Namespace Reference

Classes

- class [MSwitch](#)

5.4.1 Detailed Description

Fluidic Handling Software
Ashutosh Agarwal Lab
University of Miami

by:
Liev Birman
Adiel Hernandez

5.5 Pump Namespace Reference

Classes

- class [ThreePump](#)

Functions

- def [get_ports](#)()

5.5.1 Detailed Description

Fluidic Handling Software
Ashutosh Agarwal Lab
University of Miami

by:
Liev Birman
Adiel Hernandez

5.5.2 Function Documentation

5.5.2.1 `get_ports()`

```
def Pump.get_ports ( )
```

This method uses serial module's `comports` command to get and return serial ports.

5.6 TwoSwitch Namespace Reference

Classes

- class [TwoSwitch](#)

5.6.1 Detailed Description

Fluidic Handling Software
Ashutosh Agarwal Lab
University of Miami

by:
Liev Birman
Adiel Hernandez

Chapter 6

Class Documentation

6.1 gui.App Class Reference

Inherits Tk.

Public Member Functions

- def `__init__` (self, *args, **kwargs)
- def `show_frame` (self, cont)
- def `get_comport` (self, cmd, response)
- def `get_ports` (self)
- def `set_defaults` (self)
- def `done_and_load` (self)
- def `positive_run_status` (self)
- def `negative_run_status` (self)
- def `message_window` (self, message)

Public Attributes

- `appHighlightFont`
- `buttonFont`
- `pageFont`
- `ports`
- `max_flowrate_1`
- `max_flowrate_2`
- `max_flowrate_3`
- `max_flowrate_list`
- `default_diameter_index`
- `well_working_volume`
- `hasPump`
- *DEVICE HANDLING #####.*
- `hasMani`
- `hasColl`
- `portPump`
- `portMani`
- `portColl`

- [port2Switch](#)
- [myPump](#)
- [myMani](#)
- [myColl](#)
- [my2Switch](#)
- [should_be_running](#)
- [frames](#)

FRAME HANDLING #####.

6.1.1 Detailed Description

The App class

- Does window handling
- Holds objects and variables that need to be referenced by multiple windows

6.1.2 Constructor & Destructor Documentation

6.1.2.1 `__init__()`

```
def gui.App.__init__ (
    self,
    * args,
    ** kwargs )
```

- define window size and expandability
- define fonts
- instantiates windows

6.1.3 Member Function Documentation

6.1.3.1 `done_and_load()`

```
def gui.App.done_and_load (
    self )
```

This function checks if we have these devices and then creates instances for each one. Open ManualPage when done. When we've finished specifying our devices in the WelcomePage window we run this by clicking done.

6.1.3.2 get_comport()

```
def gui.App.get_comport (
    self,
    cmd,
    response )
```

Use this function to automatically identify a device. Input a command the device recognizes and a response you know it should return. Make sure to use unique commands when spotting multiple devices.

6.1.3.3 get_ports()

```
def gui.App.get_ports (
    self )
```

fills self.ports list with comports currently connected to device

6.1.3.4 message_window()

```
def gui.App.message_window (
    self,
    message )
```

Creates a pop-up message of your choice formatted with formatting from App class

6.1.3.5 negative_run_status()

```
def gui.App.negative_run_status (
    self )
```

sets self.should_be_running to False.

6.1.3.6 positive_run_status()

```
def gui.App.positive_run_status (
    self )

sets self.should_be_running to True.
```

6.1.3.7 set_defaults()

```
def gui.App.set_defaults (
    self )

Sets defaults as defined within function
```

6.1.3.8 show_frame()

```
def gui.App.show_frame (
    self,
    cont )

Raises frame from self.frames dictionary
```

6.1.4 Member Data Documentation

6.1.4.1 appHighlightFont

gui.App.appHighlightFont

6.1.4.2 buttonFont

gui.App.buttonFont

6.1.4.3 default_diameter_index

`gui.App.default_diameter_index`

6.1.4.4 frames

`gui.App.frames`

FRAME HANDLING #####.

6.1.4.5 hasColl

`gui.App.hasColl`

6.1.4.6 hasMani

`gui.App.hasMani`

6.1.4.7 hasPump

`gui.App.hasPump`

DEVICE HANDLING #####.

6.1.4.8 max_flowrate_1

`gui.App.max_flowrate_1`

6.1.4.9 max_flowrate_2

`gui.App.max_flowrate_2`

6.1.4.10 max_flowrate_3

`gui.App.max_flowrate_3`

6.1.4.11 max_flowrate_list

`gui.App.max_flowrate_list`

6.1.4.12 my2Switch

`gui.App.my2Switch`

6.1.4.13 myColl

`gui.App.myColl`

6.1.4.14 myMani

`gui.App.myMani`

6.1.4.15 myPump

`gui.App.myPump`

6.1.4.16 pageFont

`gui.App.pageFont`

6.1.4.17 port2Switch

`gui.App.port2Switch`

6.1.4.18 portColl

`gui.App.portColl`

6.1.4.19 portMani

`gui.App.portMani`

6.1.4.20 portPump

`gui.App.portPump`

6.1.4.21 ports

`gui.App.ports`

6.1.4.22 should_be_running

`gui.App.should_be_running`

6.1.4.23 well_working_volume

`gui.App.well_working_volume`

The documentation for this class was generated from the following file:

- [gui.py](#)

6.2 gui.**AutomaticPage** Class Reference

Inherits `Frame`.

Public Member Functions

- `def __init__ (self, parent, controller)`
- `def file_save (self)`
- `def file_load (self)`
- `def delete_step (self, step_number)`
- `def clear_recipe (self)`
- `def add_step (self)`
- `def create_step (self, channels, reservoirs)`
- `def check_channel (self, controller, channel, frame_object)`
- `def run_recipe (self, controller)`

Public Attributes

- `mainlabel`
- `channels`
- `reservoirs`
- `steplist`
- `stepnumbers`
- `numberlist`
- `number_deleted`
- `todelete`
- `deletestepcombo`
- `newstep`
- `deletestep`
- `manualbutton`
- `settingsbutton`
- `savebutton`
- `ejectbutton`
- `loadbutton`
- `runbutton`
- `clearbutton`
- `headinglabel`
- `twoSwitchValues`

6.2.1 Constructor & Destructor Documentation

6.2.1.1 __init__()

```
def gui.AutomaticPage.__init__ (
    self,
    parent,
    controller )
```

6.2.2 Member Function Documentation

6.2.2.1 add_step()

```
def gui.AutomaticPage.add_step (
    self )
```

6.2.2.2 check_channel()

```
def gui.AutomaticPage.check_channel (
    self,
    controller,
    channel,
    frame_object )
```

6.2.2.3 clear_recipe()

```
def gui.AutomaticPage.clear_recipe (
    self )
```

6.2.2.4 create_step()

```
def gui.AutomaticPage.create_step (
    self,
    channels,
    reservoirs )
```

6.2.2.5 delete_step()

```
def gui.AutomaticPage.delete_step (
    self,
    step_number )
```

Note that step number is not an int but a tk.StringVar

6.2.2.6 file_load()

```
def gui.AutomaticPage.file_load (
    self )
```

6.2.2.7 file_save()

```
def gui.AutomaticPage.file_save (
    self )
```

6.2.2.8 run_recipe()

```
def gui.AutomaticPage.run_recipe (
    self,
    controller )
```

Build in abort functionality and progress tracker

6.2.3 Member Data Documentation

6.2.3.1 channels

`gui.AutomaticPage.channels`

6.2.3.2 clearbutton

`gui.AutomaticPage.clearbutton`

6.2.3.3 deletestep

`gui.AutomaticPage.deletestep`

6.2.3.4 deletestepcombo

`gui.AutomaticPage.deletestepcombo`

6.2.3.5 ejectbutton

`gui.AutomaticPage.ejectbutton`

6.2.3.6 headinglabel

`gui.AutomaticPage.headinglabel`

6.2.3.7 loadbutton

`gui.AutomaticPage.loadbutton`

6.2.3.8 mainlabel

`gui.AutomaticPage.mainlabel`

6.2.3.9 manualbutton

`gui.AutomaticPage.manualbutton`

6.2.3.10 newstep

`gui.AutomaticPage.newstep`

6.2.3.11 number_deleted

`gui.AutomaticPage.number_deleted`

6.2.3.12 numberlist

`gui.AutomaticPage.numberlist`

6.2.3.13 reservoirs

`gui.AutomaticPage.reservoirs`

6.2.3.14 runbutton

`gui.AutomaticPage.runbutton`

6.2.3.15 savebutton

`gui.AutomaticPage.savebutton`

6.2.3.16 settingsbutton

`gui.AutomaticPage.settingsbutton`

6.2.3.17 steplist

`gui.AutomaticPage.steplist`

6.2.3.18 stepnumbers

`gui.AutomaticPage.stepnumbers`

6.2.3.19 todelete

`gui.AutomaticPage.todelete`

6.2.3.20 twoSwitchValues

`gui.AutomaticPage.twoSwitchValues`

The documentation for this class was generated from the following file:

- [gui.py](#)

6.3 Collection.Collectador Class Reference

Public Member Functions

- `def __init__ (self, my_port)`
- `def serial_connect (self)`
- `def send (self, cmd)`
- `def chop_return (self, ret)`
- `def reset (self)`
- `def eject (self)`
- `def next_site (self)`
- `def last_site (self)`
- `def moveOneUp (self)`
- `def moveOneDown (self)`
- `def moveOneRight (self)`
- `def moveOneLeft (self)`
- `def setOrigin (self)`
- `def toggle_pattern (self)`
- `def get_info (self)`

Public Attributes

- `baud`
- `port`
- `ser`
- `serialConnected`
- `position`
- `uniqueID`
- `usingSnakePattern`
- `currentPattern`
- `connected`

6.3.1 Detailed Description

This class is used for control of stepper motors. A serial connection is established with the microcontroller

Attributes:

`baud (int)`: The baud rate the serial connection is using.
`port (string)`: the microcontrollers port for the serial connection.
`ser (serial Object)`: Instance of serial object representing the serial connection.
`serialConnected (bool)`: True/False of whether the serial connection was established.
`position (int)`: the position that the XY Stage is currently at.
`uniqueID (string)`: the unique identifier the microcontroller returns for automatic port selection.
`usingSnakePattern (bool)`: True/False whether the collection pattern is currently a snake pattern.
`currentPattern (string)`: string representation of the current pattern.

6.3.2 Constructor & Destructor Documentation

6.3.2.1 `__init__()`

```
def Collection.Collectador.__init__ (
    self,
    my_port )
```

The constructor for the Collectador class.

Parameters:

`my_port` (string): the microcontrollers port for the serial connection.

6.3.3 Member Function Documentation

6.3.3.1 `chop_return()`

```
def Collection.Collectador.chop_return (
    self,
    ret )
```

This method modifies the output of the pump which comes with a carriage return and newline at the end of itself. The carriage return and newline is cut off here.

Returns:

Output of the pump without the carriage return and newline at the end.

6.3.3.2 `eject()`

```
def Collection.Collectador.eject (
    self )
```

This method sends the command "E" to the microcontroller which is programmed to move the XY Stage to the eject position. Simultaneously also updates the current position of the XY Stage.

6.3.3.3 get_info()

```
def Collection.Collectador.get_info (
    self )
```

This method sends the command "?" to the microcontroller which is programmed to send back a unique ID. The attribute uniqueID is set to the microcontrollers response.

Returns:

response (string): The unique ID that the microcontroller sends back through the serial connection.

6.3.3.4 last_site()

```
def Collection.Collectador.last_site (
    self )
```

This method sends the command "L" to the microcontroller which is programmed to move the XY Stage to the last. Simultaneously also updates the current position of the XY Stage.

6.3.3.5 moveOneDown()

```
def Collection.Collectador.moveOneDown (
    self )
```

This method makes the top motor move one step down.

6.3.3.6 moveOneLeft()

```
def Collection.Collectador.moveOneLeft (
    self )
```

This method makes the bottom motor move one step to the left.

6.3.3.7 moveOneRight()

```
def Collection.Collectador.moveOneRight (
    self )
```

This method makes the bottom motor move one step to the right.

6.3.3.8 moveOneUp()

```
def Collection.Collectador.moveOneUp (
    self )
```

This method makes the top motor move one step up.

6.3.3.9 next_site()

```
def Collection.Collectador.next_site (
    self )
```

This method sends the command "N" to the microcontroller which is programmed to move the XY Stage to the next. Simultaneously also updates the current position of the XY Stage.

6.3.3.10 reset()

```
def Collection.Collectador.reset (
    self )
```

This method sends the command "Z" to the microcontroller which is programmed to move the XY Stage to the origin. Simultaneously also updates the current position of the XY Stage.

6.3.3.11 send()

```
def Collection.Collectador.send (
    self,
    cmd )
```

This method sends a command across the serial connection.

Parameters:

cmd (string): The command or string that is to be sent to the microcontroller.

6.3.3.12 serial_connect()

```
def Collection.Collectador.serial_connect (
    self )
```

This method establishes the serial connection with the microcontroller.

Once the comport of the pump is known we open a serial connection to it using pySerial.

6.3.3.13 setOrigin()

```
def Collection.Collectador.setOrigin (
    self )
```

This method sets the origin once the user has manually moved it to his liking.
Should be the upper left corner of well plate.

6.3.3.14 toggle_pattern()

```
def Collection.Collectador.toggle_pattern (
    self )
```

This method sends the command "X" to the microcontroller which is programmed to toggle the movement layout.
The two movement patterns currently supported are snake pattern and top down pattern.

6.3.4 Member Data Documentation

6.3.4.1 baud

Collection.Collectador.baud

6.3.4.2 connected

Collection.Collectador.connected

6.3.4.3 `currentPattern`

`Collection.Collectador.currentPattern`

6.3.4.4 `port`

`Collection.Collectador.port`

6.3.4.5 `position`

`Collection.Collectador.position`

6.3.4.6 `ser`

`Collection.Collectador.ser`

6.3.4.7 `serialConnected`

`Collection.Collectador.serialConnected`

6.3.4.8 `uniqueID`

`Collection.Collectador.uniqueID`

6.3.4.9 `usingSnakePattern`

`Collection.Collectador.usingSnakePattern`

The documentation for this class was generated from the following file:

- [Collection.py](#)

6.4 gui.ManualPage Class Reference

Inherits Frame.

Public Member Functions

- def `__init__` (self, parent, controller)
- def `create_channels` (self, controller, `channels`, `reservoirs`)
- def `stop` (self, controller, channel)
- def `stop_all` (self, controller)
- def `set_and_start` (self, controller, channel)
- def `set_and_start_all` (self, controller)
- def `start_and_collect` (self, controller, `sampling_rate`)
- def `next` (self, controller)
- def `prev` (self, controller)
- def `reset` (self, controller)
- def `eject` (self, controller)
- def `toggle` (self, controller)
- def `setOrigin` (self, controller)
- def `moveOneUp` (self, controller)
- def `moveOneDown` (self, controller)
- def `moveOneRight` (self, controller)
- def `moveOneLeft` (self, controller)
- def `done` (self, controller, window)
- def `check_all_manual_and_sampling` (self, controller)
- def `check_all_manual` (self, controller)
- def `check_channel` (self, controller, channel)

Public Attributes

- `channels`
- `reservoirs`
- `samplingrate`
- `sites`
- `autopagebutton`
- `settingsbutton`
- `startallbutton`
- `startandcollectbutton`
- `stopallbutton`
- `nextwellbutton`
- `lastwellbutton`
- `resetbutton`
- `ejectbutton`
- `togglebutton`
- `setOriginButton`
- `samplingrateentry`
- `sitesentry`
- `mainlabel`
- `autocollectlabel`
- `reservoirlabel`
- `sampleslabel`
- `CollectionControlLabel`

- [positionLabel](#)
- [patternLabel](#)
- [reservoircombobox](#)
- [greaterlabelframe](#)
- [channeltextlist](#)
- [labelframelist](#)
- [flowrateentrylist](#)
- [startbuttonlist](#)
- [stopbuttonlist](#)
- [reservoirvalues](#)
- [toplabels](#)
- [twoSwitchComboBoxList](#)
- [twoSwitchValues](#)

6.4.1 Constructor & Destructor Documentation

6.4.1.1 `__init__()`

```
def gui.ManualPage.__init__ (
    self,
    parent,
    controller )
```

6.4.2 Member Function Documentation

6.4.2.1 `check_all_manual()`

```
def gui.ManualPage.check_all_manual (
    self,
    controller )
```

use when running all channels together

6.4.2.2 `check_all_manual_and_sampling()`

```
def gui.ManualPage.check_all_manual_and_sampling (
    self,
    controller )
```

use when running all channels together

6.4.2.3 check_channel()

```
def gui.ManualPage.check_channel (
    self,
    controller,
    channel )
```

6.4.2.4 create_channels()

```
def gui.ManualPage.create_channels (
    self,
    controller,
    channels,
    reservoirs )
```

creates each channel for manual page. properties of channels are stored in lists in self. channels and reservoirs variables refer to res and ch number

6.4.2.5 done()

```
def gui.ManualPage.done (
    self,
    controller,
    window )
```

6.4.2.6 eject()

```
def gui.ManualPage.eject (
    self,
    controller )
```

6.4.2.7 moveOneDown()

```
def gui.ManualPage.moveOneDown (
    self,
    controller )
```

6.4.2.8 moveOneLeft()

```
def gui.ManualPage.moveOneLeft (
    self,
    controller )
```

6.4.2.9 moveOneRight()

```
def gui.ManualPage.moveOneRight (
    self,
    controller )
```

6.4.2.10 moveOneUp()

```
def gui.ManualPage.moveOneUp (
    self,
    controller )
```

6.4.2.11 next()

```
def gui.ManualPage.next (
    self,
    controller )
```

6.4.2.12 prev()

```
def gui.ManualPage.prev (
    self,
    controller )
```

6.4.2.13 reset()

```
def gui.ManualPage.reset (
    self,
    controller )
```


6.4.2.14 set_and_start()

```
def gui.ManualPage.set_and_start (
    self,
    controller,
    channel )
```

sends input direction, flowrate, and reservoir to pump
starts specified pump channel

6.4.2.15 set_and_start_all()

```
def gui.ManualPage.set_and_start_all (
    self,
    controller )
```

sets all of the channels and then starts pump

6.4.2.16 setOrigin()

```
def gui.ManualPage.setOrigin (
    self,
    controller )
```

6.4.2.17 start_and_collect()

```
def gui.ManualPage.start_and_collect (
    self,
    controller,
    sampling_rate )
```

6.4.2.18 stop()

```
def gui.ManualPage.stop (
    self,
    controller,
    channel )
```

stops pump channel

6.4.2.19 stop_all()

```
def gui.ManualPage.stop_all (  
    self,  
    controller )
```

```
stops pump channel
```

6.4.2.20 toggle()

```
def gui.ManualPage.toggle (  
    self,  
    controller )
```

6.4.3 Member Data Documentation

6.4.3.1 autocollectlabel

```
gui.ManualPage.autocollectlabel
```

6.4.3.2 autopagebutton

```
gui.ManualPage.autopagebutton
```

6.4.3.3 channels

```
gui.ManualPage.channels
```

6.4.3.4 channeltextlist

```
gui.ManualPage.channeltextlist
```

6.4.3.5 CollectionControlLabel

`gui.ManualPage.CollectionControlLabel`

6.4.3.6 ejectbutton

`gui.ManualPage.ejectbutton`

6.4.3.7 flowrateentrylist

`gui.ManualPage.flowrateentrylist`

6.4.3.8 greaterlabelframe

`gui.ManualPage.greaterlabelframe`

6.4.3.9 labelframelist

`gui.ManualPage.labelframelist`

6.4.3.10 lastwellbutton

`gui.ManualPage.lastwellbutton`

6.4.3.11 mainlabel

`gui.ManualPage.mainlabel`

6.4.3.12 nextwellbutton

`gui.ManualPage.nextwellbutton`

6.4.3.13 patternLabel

`gui.ManualPage.patternLabel`

6.4.3.14 positionLabel

`gui.ManualPage.positionLabel`

6.4.3.15 reservoircombobox

`gui.ManualPage.reservoircombobox`

6.4.3.16 reservoirlabel

`gui.ManualPage.reservoirlabel`

6.4.3.17 reservoirs

`gui.ManualPage.reservoirs`

6.4.3.18 reservoirvalues

`gui.ManualPage.reservoirvalues`

6.4.3.19 resetbutton

`gui.ManualPage.resetbutton`

6.4.3.20 sampleslabel

`gui.ManualPage.sampleslabel`

6.4.3.21 samplingrate

`gui.ManualPage.samplingrate`

6.4.3.22 samplingrateentry

`gui.ManualPage.samplingrateentry`

6.4.3.23 setOriginButton

`gui.ManualPage.setOriginButton`

6.4.3.24 settingsbutton

`gui.ManualPage.settingsbutton`

6.4.3.25 sites

`gui.ManualPage.sites`

6.4.3.26 sitesentry

`gui.ManualPage.sitesentry`

6.4.3.27 startallbutton

`gui.ManualPage.startallbutton`

6.4.3.28 startandcollectbutton

`gui.ManualPage.startandcollectbutton`

6.4.3.29 startbuttonlist

`gui.ManualPage.startbuttonlist`

6.4.3.30 stopallbutton

`gui.ManualPage.stopallbutton`

6.4.3.31 stopbuttonlist

`gui.ManualPage.stopbuttonlist`

6.4.3.32 togglebutton

`gui.ManualPage.togglebutton`

6.4.3.33 toplabels

`gui.ManualPage.toplabels`

6.4.3.34 twoSwitchComboBoxList

`gui.ManualPage.twoSwitchComboBoxList`

6.4.3.35 twoSwitchValues

`gui.ManualPage.twoSwitchValues`

The documentation for this class was generated from the following file:

- [gui.py](#)

6.5 Mswitch.MSwitch Class Reference

Public Member Functions

- def `__init__` (self, my_port)
- def `serial_connect` (self)
- def `send` (self, cmd)
- def `chop_return` (self, ret)
- def `set_reservoir` (self, res)
- def `get_info` (self)

Public Attributes

- `baud`
- `port`
- `ser`
- `serialConnected`
- `res`
- `uniqueID`
- `connected`

6.5.1 Detailed Description

This class is used for control of the Manifold Switch. A serial connection is established with the microcontroller.

Attributes:

```
baud (int): The baud rate the serial connection is using.  
port (string): the microcontrollers port for the serial connection.  
ser (serial Object): Instance of serial object representing the serial connection.  
serialConnected (bool): True/False of whether the serial connection was established.  
res (int): The reservoir chosen.  
uniqueID (string): the unique identifier the microcontroller returns for automatic port selection.
```

6.5.2 Constructor & Destructor Documentation

6.5.2.1 `__init__()`

```
def Mswitch.MSwitch.__init__ (  
    self,  
    my_port )
```

The constructor for the MSwitch class.

Parameters:

```
my_port (string): the microcontrollers port for the serial connection.
```

6.5.3 Member Function Documentation

6.5.3.1 chop_return()

```
def Mswitch.MSwitch.chop_return (
    self,
    ret )
```

This method modifies the output of the pump which comes with a carriage return and newline at the end of itself. The carriage return and newline is cut off here.

Returns:

Output of the pump without the carriage return and newline at the end.

6.5.3.2 get_info()

```
def Mswitch.MSwitch.get_info (
    self )
```

This method sends the command "?" to the microcontroller which is programmed to send back a unique ID. The attribute

Returns:

response (string): The unique ID that the microcontroller sends back through the serial connection.

6.5.3.3 send()

```
def Mswitch.MSwitch.send (
    self,
    cmd )
```

This method sends a command across the serial connection.

Parameters:

cmd (string): The command or string that is to be sent to the microcontroller.

6.5.3.4 serial_connect()

```
def Mswitch.MSwitch.serial_connect (
    self )
```

This method establishes the serial connection with the microcontroller.

Once the comport of the pump is known we open a serial connection to it using pySerial.

6.5.3.5 set_reservoir()

```
def Mswitch.MSwitch.set_reservoir (
    self,
    res )
```

This method sets the reservoir to pull fluid from.

6.5.4 Member Data Documentation

6.5.4.1 baud

Mswitch.MSwitch.baud

6.5.4.2 connected

Mswitch.MSwitch.connected

6.5.4.3 port

Mswitch.MSwitch.port

6.5.4.4 res

Mswitch.MSwitch.res

6.5.4.5 ser

Mswitch.MSwitch.ser

6.5.4.6 serialConnected

`Mswitch.MSwitch.serialConnected`

6.5.4.7 uniqueID

`Mswitch.MSwitch.uniqueID`

The documentation for this class was generated from the following file:

- [Mswitch.py](#)

6.6 fakeSerial.Serial Class Reference

Public Member Functions

- `def __init__ (self, port='COM1', baudrate=19200, timeout=1, bytesize=8, parity='N', stopbits=1, xonxoff=0, rtscts=0)`
- `def isOpen (self)`
- `def open (self)`
- `def close (self)`
- `def write (self, string)`
- `def read (self, n=1)`
- `def flush (self)`
- `def flushInput (self)`
- `def flushOutput (self)`
- `def readline (self)`
- `def __str__ (self)`

Public Attributes

- `name`
- `port`
- `timeout`
- `parity`
- `baudrate`
- `bytesize`
- `stopbits`
- `xonxoff`
- `rtscts`

6.6.1 Detailed Description

This class is used for simulation of an Arduino Serial Connection.

Attributes:

- name (string): the name of the port.
- port (string): the microcontrollers port for the serial connection.
- baudrate (int): The baud rate the serial connection is using.
- timeout (int): maximum milliseconds to wait for serial data.
- parity (string): parity bit type.
- stopbits (int): amount of stop bits.
- xonxoff (int): whether xonxoff is on or off for flow control.
- rtscts (int): whether rtscts flow control is on or off.

6.6.2 Constructor & Destructor Documentation

6.6.2.1 __init__()

```
def fakeSerial.Serial.__init__ (
    self,
    port = 'COM1',
    baudrate = 19200,
    timeout = 1,
    bytesize = 8,
    parity = 'N',
    stopbits = 1,
    xonxoff = 0,
    rtscts = 0 )
```

The constructor for the Serial class.

Parameters:

- port (string): the microcontrollers port for the serial connection. Defaulted to 'COM1'.
- baudrate (int): The baud rate the serial connection is using. Defaulted to 19200.
- timeout (int): maximum milliseconds to wait for serial data. Defaulted to 1.
- bytesize (int): The size of bytes. Defaulted to 8.
- parity (string): parity bit type. Defaulted to 'N'.
- stopbits (int): amount of stop bits. Defaulted to 1.
- xonxoff (int): whether xonxoff is on or off for flow control. Defaulted to 0.
- rtscts (int): whether rtscts flow control is on or off. Defaulted to 0

6.6.3 Member Function Documentation

6.6.3.1 __str__()

```
def fakeSerial.Serial.__str__ (
    self )
```

This method creates a string representation of the Serial class.

Returns:

- string: String representation of Serial class.

6.6.3.2 close()

```
def fakeSerial.Serial.close (  
    self )
```

This method closes the port artificially.

6.6.3.3 flush()

```
def fakeSerial.Serial.flush (  
    self )
```

This method flushes the port artificially.

6.6.3.4 flushInput()

```
def fakeSerial.Serial.flushInput (  
    self )
```

This method flushes the input of the port artificially.

6.6.3.5 flushOutput()

```
def fakeSerial.Serial.flushOutput (  
    self )
```

This method flushes the output of the port artificially.

6.6.3.6 isOpen()

```
def fakeSerial.Serial.isOpen (  
    self )
```

This method returns True if the port to the arduino is open. False otherwise.

Returns:

 _isOpen (bool): True/False depending if serial connection is open.

6.6.3.7 open()

```
def fakeSerial.Serial.open (
    self )
```

This method opens the port artificially.

6.6.3.8 read()

```
def fakeSerial.Serial.read (
    self,
    n = 1 )
```

This method simulates reading a pretermind amount of characters from the Arduino.

Parameters:

n (int): number of characters to read.

Returns:

s (string): the characters read from the Arduino.

6.6.3.9 readline()

```
def fakeSerial.Serial.readline (
    self )
```

This method simulates reading an entire line of characters until a "\n" is encountered.

Returns:

s (string): The line of characters read if there is one. Otherwise, an empty string.

6.6.3.10 write()

```
def fakeSerial.Serial.write (
    self,
    string )
```

This method simulates writing a string of characters to the Arduino.

Returns:

An empty string.

6.6.4 Member Data Documentation

6.6.4.1 baudrate

`fakeSerial.Serial.baudrate`

6.6.4.2 bytesize

`fakeSerial.Serial.bytesize`

6.6.4.3 name

`fakeSerial.Serial.name`

6.6.4.4 parity

`fakeSerial.Serial.parity`

6.6.4.5 port

`fakeSerial.Serial.port`

6.6.4.6 rtscts

`fakeSerial.Serial.rtscts`

6.6.4.7 stopbits

`fakeSerial.Serial.stopbits`

6.6.4.8 timeout

```
fakeSerial.Serial.timeout
```

6.6.4.9 xonxoff

```
fakeSerial.Serial.xonxoff
```

The documentation for this class was generated from the following file:

- [fakeSerial.py](#)

6.7 gui.SettingsPage Class Reference

Inherits Frame.

Public Member Functions

- def [__init__](#) (self, parent, controller)
- def [update](#) (self, controller)
- def [calibrate](#) (self, controller)

Public Attributes

- [channels](#)
- [reservoirs](#)
- [mainlabel](#)
- [autopagebutton](#)
- [manualbutton](#)
- [updatebutton](#)
- [calibratecombobox](#)
- [channelvalues](#)
- [calibratebutton](#)
- [measuredVolume](#)
- [measuredVolumeEntry](#)
- [entercalibratedbutton](#)
- [calibrationlabel](#)
- [channelcaliblabel](#)
- [directionheading](#)
- [directionslist](#)
- [directionlabellist](#)
- [directions](#)
- [diameterlist](#)
- [diameterlabel](#)
- [diametervalues](#)

6.7.1 Detailed Description

Everything that is being sent and all experimental parameters should be recorded here. Start and stop events. Serial commands. Status.

6.7.2 Constructor & Destructor Documentation

6.7.2.1 `__init__()`

```
def gui.SettingsPage.__init__ (
    self,
    parent,
    controller )
```

6.7.3 Member Function Documentation

6.7.3.1 `calibrate()`

```
def gui.SettingsPage.calibrate (
    self,
    controller )
```

6.7.3.2 `update()`

```
def gui.SettingsPage.update (
    self,
    controller )
```

6.7.4 Member Data Documentation

6.7.4.1 `autopagebutton`

```
gui.SettingsPage.autopagebutton
```


6.7.4.2 calibratebutton

`gui.SettingsPage.calibratebutton`

6.7.4.3 calibratecombobox

`gui.SettingsPage.calibratecombobox`

6.7.4.4 calibrationlabel

`gui.SettingsPage.calibrationlabel`

6.7.4.5 channelcaliblabel

`gui.SettingsPage.channelcaliblabel`

6.7.4.6 channels

`gui.SettingsPage.channels`

6.7.4.7 channelvalues

`gui.SettingsPage.channelvalues`

6.7.4.8 diameterlabel

`gui.SettingsPage.diameterlabel`

6.7.4.9 diameterlist

`gui.SettingsPage.diameterlist`

6.7.4.10 diametervalues

`gui.SettingsPage.diametervalues`

6.7.4.11 directionheading

`gui.SettingsPage.directionheading`

6.7.4.12 directionlabellist

`gui.SettingsPage.directionlabellist`

6.7.4.13 directions

`gui.SettingsPage.directions`

6.7.4.14 directionslist

`gui.SettingsPage.directionslist`

6.7.4.15 entercalibratedbutton

`gui.SettingsPage.entercalibratedbutton`

6.7.4.16 mainlabel

`gui.SettingsPage.mainlabel`

6.7.4.17 manualbutton

`gui.SettingsPage.manualbutton`

6.7.4.18 measuredVolume

`gui.SettingsPage.measuredVolume`

6.7.4.19 measuredVolumeEntry

`gui.SettingsPage.measuredVolumeEntry`

6.7.4.20 reservoirs

`gui.SettingsPage.reservoirs`

6.7.4.21 updatebutton

`gui.SettingsPage.updatebutton`

The documentation for this class was generated from the following file:

- [gui.py](#)

6.8 Pump.ThreePump Class Reference

Public Member Functions

- `def __init__ (self, my_port)`
- `def serial_connect (self)`
- `def send (self, cmd)`
- `def send_return (self, cmd)`
- `def chop_return (self, ret)`
- `def FormatVolume (self, V, unit)`
- `def setFlow (self, channel, flowrate)`
- `def setDir (self, channel, direction)`
- `def start (self, channel)`
- `def stop (self, channel)`
- `def calibrate (self, channel)`
- `def abort_calibration (self, channel)`
- `def setTargetCalibrationVolume (self)`
- `def setMeasuredVolume (self, channel, volume)`
- `def start_all (self)`
- `def stop_all (self)`
- `def setDefaults (self)`

Public Attributes

- [baud](#)
- [port](#)
- [ser](#)
- [serialConnected](#)
- [calibrationvolume](#)
- [calibrationunit](#)
- [isOn](#)
- [connected](#)

6.8.1 Detailed Description

This class is used for control of the Pump. A serial connection is established with the microcontroller control

Attributes:

```
baud (int): The baud rate the serial connection is using.
port (string): the microcontrollers port for the serial connection.
ser (serial Object): Instance of serial object representing the serial connection.
serialConnected (bool): True/False of whether the serial connection was established.
calibrationvolume (int): amount of fluid used for calibration.
calibrationunit (string): the unit for the amount of fluid used for calibration.
isOn (bool): True/False on whether the pump is on or not.
```

6.8.2 Constructor & Destructor Documentation

6.8.2.1 __init__()

```
def Pump.ThreePump.__init__ (
    self,
    my_port )
```

The constructor for the TwoSwitch class.

Parameters:

```
my_port (string): the microcontrollers port for the serial connection.
```

6.8.3 Member Function Documentation

6.8.3.1 abort_calibration()

```
def Pump.ThreePump.abort_calibration (
    self,
    channel )
```

This method sends a command to the microcontroller which is programmed to abort calibration of the pump only f

Parameters:

```
channel (int): the channel that will abort calibration.
```

6.8.3.2 `calibrate()`

```
def Pump.ThreePump.calibrate (
    self,
    channel )
```

This method sends a command to the microcontroller which is programmed to calibrate the pump only for the channel.

Parameters:

channel (int): the channel that will be calibrated.

6.8.3.3 `chop_return()`

```
def Pump.ThreePump.chop_return (
    self,
    ret )
```

This method modifies the output of the two switch which comes with a carriage return and newline at the end of the line. The carriage return and newline is cut off here.

Returns:

Output of the pump without the carriage return and newline at the end.

6.8.3.4 `FormatVolume()`

```
def Pump.ThreePump.FormatVolume (
    self,
    V,
    unit )
```

This method formats the volume.

Parameters:

V (int): The amount of volume.

unit (string): The unit the volume should be set in.

6.8.3.5 `send()`

```
def Pump.ThreePump.send (
    self,
    cmd )
```

This method sends a command across the serial connection.

Parameters:

cmd (string): The command or string that is to be sent to the microcontroller.

6.8.3.6 send_return()

```
def Pump.ThreePump.send_return (
    self,
    cmd )
```

This method uses the serial connection opened instance of pump and sends the text written in cmd across that connection.

Returns:

response (string): The response or what the microcontroller returns.

6.8.3.7 serial_connect()

```
def Pump.ThreePump.serial_connect (
    self )
```

This method establishes the serial connection with the microcontroller.

Once the comport of the pump is known we open a serial connection to it using pySerial.

6.8.3.8 setDefaults()

```
def Pump.ThreePump.setDefaults (
    self )
```

This method sets the defaults for the pump.

6.8.3.9 setDir()

```
def Pump.ThreePump.setDir (
    self,
    channel,
    direction )
```

This method sets the direction for the flow of fluid.

Parameters:

channel (int): The channel flow rate direction will be set.

direction (string): The direction for the pump CW or CCW.

6.8.3.10 setFlow()

```
def Pump.ThreePump.setFlow (
    self,
    channel,
    flowrate )
```

This method sets the flow rate of the pump.

Parameters:

channel (int): The channel whose flow rate will be set.
flowrate (float): The flow rate which is to be set.

6.8.3.11 setMeasuredVolume()

```
def Pump.ThreePump.setMeasuredVolume (
    self,
    channel,
    volume )
```

This method sets the measured volume.

6.8.3.12 setTargetCalibrationVolume()

```
def Pump.ThreePump.setTargetCalibrationVolume (
    self )
```

This method sets the target calibration volume.

6.8.3.13 start()

```
def Pump.ThreePump.start (
    self,
    channel )
```

This method sends a command to the microcontroller which is programmed to start the pump only for the channel

Parameters:

channel (int): the channel that will start.

6.8.3.14 start_all()

```
def Pump.ThreePump.start_all (
    self )
```

This method sends a command to the microcontroller which is programmed to start the pump for all channels.

6.8.3.15 stop()

```
def Pump.ThreePump.stop (
    self,
    channel )
```

This method sends a command to the microcontroller which is programmed to stop the pump only for the channel *channel*.

Parameters:

channel (int): the channel that will stop.

6.8.3.16 stop_all()

```
def Pump.ThreePump.stop_all (
    self )
```

This method sends a command to the microcontroller which is programmed to stop the pump for all the channels.

6.8.4 Member Data Documentation

6.8.4.1 baud

Pump.ThreePump.baud

6.8.4.2 calibrationunit

Pump.ThreePump.calibrationunit

6.8.4.3 calibrationvolume

`Pump.ThreePump.calibrationvolume`

6.8.4.4 connected

`Pump.ThreePump.connected`

6.8.4.5 isOn

`Pump.ThreePump.isOn`

6.8.4.6 port

`Pump.ThreePump.port`

6.8.4.7 ser

`Pump.ThreePump.ser`

6.8.4.8 serialConnected

`Pump.ThreePump.serialConnected`

The documentation for this class was generated from the following file:

- [Pump.py](#)

6.9 TwoSwitch.TwoSwitch Class Reference

Public Member Functions

- def [__init__](#) (self, my_port)
- def [serial_connect](#) (self)
- def [send](#) (self, cmd)
- def [chop_return](#) (self, ret)
- def [setRecirculate](#) (self, channel)
- def [setCollect](#) (self, channel)
- def [getIdentifier](#) (self)

Public Attributes

- [baud](#)
- [port](#)
- [ser](#)
- [serialConnected](#)
- [uniqueID](#)
- [willRecirculate](#)
- [connected](#)

6.9.1 Detailed Description

This class is used for control of the two switches. A serial connection is established with the microcontroller.

Attributes:

`baud (int)`: The baud rate the serial connection is using.
`port (string)`: the microcontrollers port for the serial connection.
`ser (serial Object)`: Instance of serial object representing the serial connection.
`serialConnected (bool)`: True/False of whether the serial connection was established.
`uniqueID (string)`: the unique identifier the microcontroller returns for automatic port selection.
`willRecirculate (bool)`: True/False as to if the two switches are set to recirculate.

6.9.2 Constructor & Destructor Documentation

6.9.2.1 `__init__()`

```
def TwoSwitch.TwoSwitch.__init__ (
    self,
    my_port )
```

The constructor for the TwoSwitch class.

Parameters:

`my_port (string)`: the microcontrollers port for the serial connection.

6.9.3 Member Function Documentation

6.9.3.1 `chop_return()`

```
def TwoSwitch.TwoSwitch.chop_return (
    self,
    ret )
```

This method modifies the output of the two switch which comes with a carriage return and newline at the end of the output. The carriage return and newline is cut off here.

Returns:

Output of the pump without the carriage return and newline at the end.

6.9.3.2 getIdentifier()

```
def TwoSwitch.TwoSwitch.getIdentifier (
    self )
```

This method sends the command "?" to the microcontroller which is programmed to send back a unique ID. The att

Returns:

response (string): The unique ID that the microcontroller sends back through the serial connection.

6.9.3.3 send()

```
def TwoSwitch.TwoSwitch.send (
    self,
    cmd )
```

This method sends a command across the serial connection.

Parameters:

cmd (string): The command or string that is to be sent to the microcontroller.

6.9.3.4 serial_connect()

```
def TwoSwitch.TwoSwitch.serial_connect (
    self )
```

This method establishes the serial connection with the microcontroller.

Once the comport of the pump is known we open a serial connection to it using pySerial.

6.9.3.5 setCollect()

```
def TwoSwitch.TwoSwitch.setCollect (
    self,
    channel )
```

This method sends a command to the microcontroller based on what channel needs to be set to recirculate.

Parameters:

channel (int): The fluidic channel chosen to be set to collect.

6.9.3.6 setRecirculate()

```
def TwoSwitch.TwoSwitch.setRecirculate (
    self,
    channel )
```

This method sends a command to the microcontroller based on what channel needs to be set to recirculate.

Parameters:

channel (int): The fluidic channel chosen to be set to recirculate.

6.9.4 Member Data Documentation

6.9.4.1 baud

TwoSwitch.TwoSwitch.baud

6.9.4.2 connected

TwoSwitch.TwoSwitch.connected

6.9.4.3 port

TwoSwitch.TwoSwitch.port

6.9.4.4 ser

TwoSwitch.TwoSwitch.ser

6.9.4.5 serialConnected

TwoSwitch.TwoSwitch.serialConnected

6.9.4.6 uniqueID

```
TwoSwitch.TwoSwitch.uniqueID
```

6.9.4.7 willRecirculate

```
TwoSwitch.TwoSwitch.willRecirculate
```

The documentation for this class was generated from the following file:

- [TwoSwitch.py](#)

6.10 gui.WelcomePage Class Reference

Inherits Frame.

Public Member Functions

- `def __init__ (self, parent, controller)`
- `def flip (self, widget)`

Public Attributes

- `mainlabel`
- `sublabelframe`
- `perfusorframe`
- `collectorframe`
- `manilabel`
- *Port & Device Selection##### finding all of the ports.*
- `pumplabel`
- `coll_label`
- `donebutton`

6.10.1 Constructor & Destructor Documentation

6.10.1.1 __init__()

```
def gui.WelcomePage.__init__ (
    self,
    parent,
    controller )
```

6.10.2 Member Function Documentation

6.10.2.1 flip()

```
def gui.WelcomePage.flip (
    self,
    widget )
```

Enables or disables a widget. It is run by the checkbuttons on the WelcomPage

6.10.3 Member Data Documentation

6.10.3.1 coll_label

```
gui.WelcomePage.coll_label
```

6.10.3.2 collectorframe

```
gui.WelcomePage.collectorframe
```

6.10.3.3 donebutton

```
gui.WelcomePage.donebutton
```

6.10.3.4 mainlabel

```
gui.WelcomePage.mainlabel
```

6.10.3.5 manilabel

```
gui.WelcomePage.manilabel
```

Port & Device Selection##### finding all of the ports.

6.10.3.6 perfusorframe

`gui.WelcomePage.perfusorframe`

6.10.3.7 pumplabel

`gui.WelcomePage.pumplabel`

6.10.3.8 sublabelframe

`gui.WelcomePage.sublabelframe`

The documentation for this class was generated from the following file:

- [gui.py](#)

Chapter 7

File Documentation

7.1 ArduinoSnakePattern.cpp File Reference

Arduino code for XYStage Movement and [TwoSwitch](#) control.

```
#include <Wire.h>
#include <Adafruit_MotorShield.h>
#include <AccelStepper.h>
#include <MultiStepper.h>
#include "utility/Adafruit_MS_PWM_ServoDriver.h"
```

Functions

- void [forwardstep1](#) ()
Wrapper Forward Step function for Stepper Motor 1.
- void [backwardstep1](#) ()
Wrapper Backward Step function for Stepper Motor 1.
- void [forwardstep2](#) ()
Wrapper Forward Step function for Stepper Motor 2.
- void [backwardstep2](#) ()
Wrapper Backward Step function for Stepper Motor 2.
- void [populateLocations](#) (long(&[locations](#))[[numLocations](#)][2])
This function populates the location matrix with the layout for moving to any location on the well plate.
- void [moveGen](#) (Adafruit_StepperMotor *motor, int delay, int myDirection, int stepType)
Generic Move Function for Stepper Motors.
- void [moveRight](#) ()
MOVE RIGHT FUNCTION.
- void [moveLeft](#) ()
MOVE LEFT FUNCTION.
- void [moveDown](#) ()
MOVE DOWN FUNCTION.
- void [moveUp](#) ()
MOVE UP FUNCTION.
- int [moveNext](#) ()
MOVE NEXT FUNCTION.

- int `moveLast` ()
MOVE NEXT FUNCTION.
- int `set` ()
- int `eject` ()
EJECT WELL PLATE FUNCTION.
- void `reset` ()
RESET FUNCTION.
- AccelStepper `stepperB` (`forwardstep1`, `backwardstep1`)
- AccelStepper `stepperT` (`forwardstep2`, `backwardstep2`)
- void `setup` ()
Arduino Setup.
- void `loop` ()
Arduino Loop Function, this runs continuously.

Variables

- const int `delayRPM` = 255
Imports needed.
- const int `rows` = 8
Amount of rows on well plate.
- const int `columns` = 12
Amount of columns on well plate.
- const int `numLocations` = ((`rows` * `columns`) / 3) + 1
Amount of triplet wells on well plate.
- const float `maximumSpeed` = 255.0
Max Speed of the Stepper Motors.
- const float `acceleration` = 2000.0
Acceleration for the Stepper Motors.
- const int `stepsPerRevolution` = 200
Steps per revolution for the Stepper Motor.
- const int `stepsPerWell` = 225
Amount of steps needed from center of well to the next.
- const int `stepsToEject` = 1000
Amount of steps to eject entire well plate.
- const int `limitSwitchT` = 37
Digital Pin Number for Top Stepper Motor Limit Switch.
- const int `limitSwitchB` = 36
Digital Pin Number for Bottom Stepper Motor Limit Switch.
- const int `TwoSwitchA` = 34
Digital Pin Number for TwoSwitchA.
- const int `TwoSwitchB` = 33
Digital Pin Number for TwoSwitchB.
- const int `TwoSwitchC` = 32
Digital Pin Number for TwoSwitchC.
- Adafruit_MotorShield `AFMS` = Adafruit_MotorShield()
- Adafruit_StepperMotor * `myMotorT` = AFMS.getStepper(`stepsPerRevolution`, 1)
- Adafruit_StepperMotor * `myMotorB` = AFMS.getStepper(`stepsPerRevolution`, 2)
- MultiStepper `steppers`
- int `myPosition` = 0
Position variable to keep track of which triplet well we are on.
- long `locations` [`numLocations`][2]
Three dimensional array for schematic of positions if MultiStepper is used.
- long `positions` [2]
Array for movement needed for using both motors in MultiStepper.

7.1.1 Detailed Description

Arduino code for XYStage Movement and [TwoSwitch](#) control.

Author

Adiel Hernandez

7.1.2 Function Documentation

7.1.2.1 backwardstep1()

```
void backwardstep1 ( )
```

Wrapper Backward Step function for Stepper Motor 1.

To be used with the AccelStepper Library.

7.1.2.2 backwardstep2()

```
void backwardstep2 ( )
```

Wrapper Backward Step function for Stepper Motor 2.

To be used with the AccelStepper Library.

7.1.2.3 eject()

```
int eject ( )
```

EJECT WELL PLATE FUNCTION.

Moves the motor accordingly so it goes to the eject well plate position. It moves the well plate down from its current location. Sets myPosition to 32.

7.1.2.4 forwardstep1()

```
void forwardstep1 ( )
```

Wrapper Forward Step function for Stepper Motor 1.

To be used with the AccelStepper Library.

7.1.2.5 forwardstep2()

```
void forwardstep2 ( )
```

Wrapper Forward Step function for Stepper Motor 2.

To be used with the AccelStepper Library.

7.1.2.6 loop()

```
void loop ( )
```

Arduino Loop Function, this runs continuously.

Essentially this looks for incoming commands which correspond to a certain function for motor movement or Limit Switching.

7.1.2.7 moveDown()

```
void moveDown ( )
```

MOVE DOWN FUNCTION.

Moves the bottom motor to the left. With SINGLE step type.

7.1.2.8 moveGen()

```
void moveGen (
    Adafruit_StepperMotor * motor,
    int delay,
    int myDirection,
    int stepType )
```

Generic Move Function for Stepper Motors.

This function takes in a pointer to the motor that has to be moved, the speed, the direction, and step type. All of this is used to move the correct motor accordingly.

Parameters

<i>*motor</i>	a Adafruit_StepperMotor object pointer for the motor to be moved.
<i>delay</i>	the speed or RPM to move the motor at.
<i>myDirection</i>	the direction to move the motor in FORWARD or BACKWARD.
<i>stepType</i>	the step type to move the motor in: SINGLE, DOUBLE, INTERLEAVE, or MICROSTEP.

7.1.2.9 moveLast()

```
int moveLast ( )
```

MOVE NEXT FUNCTION.

Moves the motor accordingly so it goes to the previous well plate given it is not on the first possible triplet well location. Moves in a snake pattern. Increments myPosition by 1.

7.1.2.10 moveLeft()

```
void moveLeft ( )
```

MOVE LEFT FUNCTION.

Moves the top motor to the left. With SINGLE step type.

7.1.2.11 moveNext()

```
int moveNext ( )
```

MOVE NEXT FUNCTION.

Moves the motor accordingly so it goes to the next well plate given it is not on the last possible triplet well location. Moves in a snake pattern. Decrements myPosition by 1.

7.1.2.12 moveRight()

```
void moveRight ( )
```

MOVE RIGHT FUNCTION.

Moves the top motor to the right. With SINGLE step type.

7.1.2.13 moveUp()

```
void moveUp ( )
```

MOVE UP FUNCTION.

Moves the bottom motor to the right. With SINGLE step type.

7.1.2.14 populateLocations()

```
void populateLocations (
    long(&) locations[numLocations][2] )
```

This function populates the location matrix with the layout for moving to any location on the well plate.

Not currently being used since we are doing a snake pattern but could be used in future.

7.1.2.15 reset()

```
void reset ( )
```

RESET FUNCTION.

Moves the motor accordingly so it goes to the origin well plate position. Moves the motors simultaneously until both limit switches are activated and then moves a predetermined amount to origin.

7.1.2.16 set()

```
int set ( )
```

7.1.2.17 setup()

```
void setup ( )
```

Arduino Setup.

This only runs once.

7.1.2.18 stepperB()

```
AccelStepper stepperB (
    forwardstep1 ,
    backwardstep1 )
```

7.1.2.19 stepperT()

```
AccelStepper stepperT (
    forwardstep2 ,
    backwardstep2 )
```

7.1.3 Variable Documentation

7.1.3.1 acceleration

```
const float acceleration = 2000.0
```

Acceleration for the Stepper Motors.

7.1.3.2 AFMS

```
Adafruit_MotorShield AFMS = Adafruit_MotorShield()
```

7.1.3.3 columns

```
const int columns = 12
```

Amount of columns on well plate.

7.1.3.4 delayRPM

```
const int delayRPM = 255
```

Imports needed.

The Adafruit_MotorShieldV2 library and AccelStepper Library is used. Global Variables. Constains mostly motor values and digital pin values. RPM for the stepper motors.

7.1.3.5 limitSwitchB

```
const int limitSwitchB = 36
```

Digital Pin Number for Bottom Stepper Motor Limit Switch.

7.1.3.6 limitSwitchT

```
const int limitSwitchT = 37
```

Digital Pin Number for Top Stepper Motor Limit Switch.

7.1.3.7 locations

```
long locations[numLocations][2]
```

Three dimensional array for schematic of positions if MultiStepper is used.

7.1.3.8 maximumSpeed

```
const float maximumSpeed = 255.0
```

Max Speed of the Stepper Motors.

7.1.3.9 myMotorB

```
Adafruit_StepperMotor* myMotorB = AFMS.getStepper(stepsPerRevolution, 2)
```

7.1.3.10 myMotorT

```
Adafruit_StepperMotor* myMotorT = AFMS.getStepper(stepsPerRevolution, 1)
```

7.1.3.11 myPosition

```
int myPosition = 0
```

Position variable to keep track of which triplet well we are on.

7.1.3.12 numLocations

```
const int numLocations = ((rows * columns) / 3) + 1
```

Amount of triplet wells on well plate.

(Plus 1 for the eject location).

7.1.3.13 positions

```
long positions[2]
```

Array for movement needed for using both motors in MultiStepper.

7.1.3.14 rows

```
const int rows = 8
```

Amount of rows on well plate.

7.1.3.15 steppers

```
MultiStepper steppers
```

7.1.3.16 stepsPerRevolution

```
const int stepsPerRevolution = 200
```

Steps per revolution for the Stepper Motor.

7.1.3.17 stepsPerWell

```
const int stepsPerWell = 225
```

Amount of steps needed from center of well to the next.

7.1.3.18 stepsToEject

```
const int stepsToEject = 1000
```

Amount of steps to eject entire well plate.

7.1.3.19 TwoSwitchA

```
const int TwoSwitchA = 34
```

Digital Pin Number for TwoSwitchA.

7.1.3.20 TwoSwitchB

```
const int TwoSwitchB = 33
```

Digital Pin Number for TwoSwitchB.

7.1.3.21 TwoSwitchC

```
const int TwoSwitchC = 32
```

Digital Pin Number for TwoSwitchC.

7.2 Collection.py File Reference

Classes

- class [Collection.Collectador](#)

Namespaces

- [Collection](#)

7.3 fakeSerial.py File Reference

Classes

- class [fakeSerial.Serial](#)

Namespaces

- [fakeSerial](#)

7.4 gui.py File Reference

Classes

- class [gui.App](#)
- class [gui.WelcomePage](#)
- class [gui.ManualPage](#)
- class [guiAutomaticPage](#)
- class [gui.SettingsPage](#)

Namespaces

- [gui](#)

Functions

- def [gui.combine_funcs](#) (*funcs)
- def [gui.message_prompt](#) (message)
- def [gui.check_input](#) (variable, desired_datatype, lower_limit=0, upper_limit=0)

Variables

- tuple [gui.LARGE_FONT](#) = ("Verdana", 10)
- [gui.app](#) = App()

7.5 locationPopulate.cpp File Reference

```
#include <iostream>
```

Functions

- void [populateLocations](#) (long(&[locations](#))[[numLocations](#)][2], int choice=1)
This function populates the location matrix with the layout for moving to any location on the well plate.
- int [main](#) ()

Variables

- const int [rows](#) = 8
- const int [columns](#) = 12
- const int [numLocations](#) = (([rows](#) * [columns](#)) / 3) + 1
- const int [stepsPerWell](#) = 225
- const int [stepsToEject](#) = 500

7.5.1 Function Documentation

7.5.1.1 main()

```
int main ( )
```

7.5.1.2 populateLocations()

```
void populateLocations (
    long(&) locations[numLocations][2],
    int choice = 1 )
```

This function populates the location matrix with the layout for moving to any location on the well plate.

Gives the user the choice for a snake pattern Choice = 1, or top to bottom pattern Choice != 1.

7.5.2 Variable Documentation

7.5.2.1 columns

```
const int columns = 12
```

7.5.2.2 numLocations

```
const int numLocations = ((rows * columns) / 3) + 1
```

7.5.2.3 rows

```
const int rows = 8
```

7.5.2.4 stepsPerWell

```
const int stepsPerWell = 225
```

7.5.2.5 stepsToEject

```
const int stepsToEject = 500
```

7.6 Mswitch.py File Reference

Classes

- class [Mswitch.MSwitch](#)

Namespaces

- [Mswitch](#)

7.7 Pump.py File Reference

Classes

- class [Pump.ThreePump](#)

Namespaces

- [Pump](#)

Functions

- def [Pump.get_ports](#) ()

7.8 TwoSwitch.py File Reference

Classes

- class [TwoSwitch.TwoSwitch](#)

Namespaces

- [TwoSwitch](#)

7.9 XYStage_LimitSwitch_TwoSwitch_Control.cpp File Reference

```
#include <Wire.h>
#include <Adafruit_MotorShield.h>
#include <AccelStepper.h>
#include <MultiStepper.h>
#include "utility/Adafruit_MS_PWM_ServoDriver.h"
```

Functions

- void `forwardstep1` ()
Wrapper Forward Step function for Stepper Motor 1.
- void `backwardstep1` ()
Wrapper Backward Step function for Stepper Motor 1.
- void `forwardstep2` ()
Wrapper Forward Step function for Stepper Motor 2.
- void `backwardstep2` ()
Wrapper Backward Step function for Stepper Motor 2.
- void `populateLocations` (long(&`locations`)[`numLocations`][2], int choice=1)
This function populates the location matrix with the layout for moving to any location on the well plate.
- void `moveGen` (Adafruit_StepperMotor *motor, int delay, int myDirection, int stepType)
Generic Move Function for Stepper Motors.
- void `moveRight` ()
MOVE RIGHT FUNCTION.
- void `moveLeft` ()
MOVE LEFT FUNCTION.
- void `moveDown` ()
MOVE DOWN FUNCTION.
- void `moveUp` ()
MOVE UP FUNCTION.
- int `moveNext` ()
MOVE NEXT FUNCTION.
- int `moveLast` ()
MOVE NEXT FUNCTION.
- int `set` ()
- int `eject` ()
EJECT WELL PLATE FUNCTION.
- void `reset` ()
RESET FUNCTION.
- AccelStepper `stepperB` (`forwardstep1`, `backwardstep1`)
- AccelStepper `stepperT` (`forwardstep2`, `backwardstep2`)
- void `setup` ()
- void `loop` ()
Arduino Loop Function, this runs continuously.

Variables

- const int `delayRPM` = 255
Imports needed.
- const int `rows` = 8
Amount of rows on well plate.
- const int `columns` = 12
Amount of columns on well plate.
- const int `numLocations` = ((`rows` * `columns`) / 3) + 1
Amount of triplet wells on well plate.
- const float `maximumSpeed` = 255.0
Max Speed of the Stepper Motors.
- const float `acceleration` = 2000.0

- Acceleration for the Stepper Motors.*
- const int `stepsPerRevolution` = 200
Steps per revolution for the Stepper Motor.
- const int `stepsPerWell` = 225
Amount of steps needed from center of well to the next.
- const int `stepsToEject` = 1000
Amount of steps to eject entire well plate.
- const int `limitSwitchT` = 41
Digital Pin Number for Top Stepper Motor Limit Switch.
- const int `limitSwitchB` = 43
Digital Pin Number for Bottom Stepper Motor Limit Switch.
- const int `TwoSwitchA` = 35
Digital Pin Number for TwoSwitchA.
- const int `TwoSwitchB` = 33
Digital Pin Number for TwoSwitchB.
- const int `TwoSwitchC` = 31
Digital Pin Number for TwoSwitchC.
- Adafruit_MotorShield `AFMS` = Adafruit_MotorShield()
- Adafruit_StepperMotor * `myMotorT` = AFMS.getStepper(`stepsPerRevolution`, 1)
- Adafruit_StepperMotor * `myMotorB` = AFMS.getStepper(`stepsPerRevolution`, 2)
- MultiStepper `steppers`
- int `myPosition` = 0
Position variable to keep track of which triplet well we are on.
- long `locations` [`numLocations`][2]
Three dimensional array for schematic of positions if MultiStepper is used.
- long `positions` [2]
Array for movement needed for using both motors in MultiStepper.
- int `collectionPattern` = 1
Collection Pattern variable to choose between snake pattern XY Stage movement (default = 1) or Top down XY Stage movement (!= 1).

7.9.1 Function Documentation

7.9.1.1 backwardstep1()

```
void backwardstep1 ( )
```

Wrapper Backward Step function for Stepper Motor 1.

To be used with the AccelStepper Library.

7.9.1.2 backwardstep2()

```
void backwardstep2 ( )
```

Wrapper Backward Step function for Stepper Motor 2.

To be used with the AccelStepper Library.

7.9.1.3 eject()

```
int eject ( )
```

EJECT WELL PLATE FUNCTION.

Moves the motor accordingly so it goes to the eject well plate position. It moves the well plate down from its current location. Sets myPosition to 32.

7.9.1.4 forwardstep1()

```
void forwardstep1 ( )
```

Wrapper Forward Step function for Stepper Motor 1.

To be used with the AccelStepper Library.

7.9.1.5 forwardstep2()

```
void forwardstep2 ( )
```

Wrapper Forward Step function for Stepper Motor 2.

To be used with the AccelStepper Library.

7.9.1.6 loop()

```
void loop ( )
```

Arduino Loop Function, this runs continuously.

Essentially this looks for incoming commands which correspond to a certain function for motor movement or Limit Switching.

7.9.1.7 moveDown()

```
void moveDown ( )
```

MOVE DOWN FUNCTION.

Moves the bottom motor to the left. With SINGLE step type.

7.9.1.8 moveGen()

```
void moveGen (
    Adafruit_StepperMotor * motor,
    int delay,
    int myDirection,
    int stepType )
```

Generic Move Function for Stepper Motors.

This function takes in a pointer to the motor that has to be moved, the speed, the direction, and step type. All of this is used to move the correct motor accordingly.

Parameters

<i>*motor</i>	a Adafruit_StepperMotor object pointer for the motor to be moved.
<i>delay</i>	the speed or RPM to move the motor at.
<i>myDirection</i>	the direction to move the motor in FORWARD or BACKWARD.
<i>stepType</i>	the step type to move the motor in: SINGLE, DOUBLE, INTERLEAVE, or MICROSTEP.

7.9.1.9 moveLast()

```
int moveLast ( )
```

MOVE NEXT FUNCTION.

Moves the motor accordingly so it goes to the previous well plate given it is not on the first possible triplet well location. Moves in a snake pattern. Increments myPosition by 1.

7.9.1.10 moveLeft()

```
void moveLeft ( )
```

MOVE LEFT FUNCTION.

Moves the top motor to the left. With SINGLE step type.

7.9.1.11 moveNext()

```
int moveNext ( )
```

MOVE NEXT FUNCTION.

Moves the motor accordingly so it goes to the next well plate given it is not on the last possible triplet well location. Moves in a snake pattern. Decrements myPosition by 1.

7.9.1.12 moveRight()

```
void moveRight ( )
```

MOVE RIGHT FUNCTION.

Moves the top motor to the right. With SINGLE step type.

7.9.1.13 moveUp()

```
void moveUp ( )
```

MOVE UP FUNCTION.

Moves the bottom motor to the right. With SINGLE step type.

7.9.1.14 populateLocations()

```
void populateLocations (
    long(&) locations[numLocations][2],
    int choice = 1 )
```

This function populates the location matrix with the layout for moving to any location on the well plate.

Gives the user the choice for a snake pattern Choice = 1, or top to bottom pattern Choice != 1.

7.9.1.15 reset()

```
void reset ( )
```

RESET FUNCTION.

Moves the motor accordingly so it goes to the origin well plate position. Moves the motors simultaneously until both limit switches are activated and then moves a predetermined amount to origin.

7.9.1.16 set()

```
int set ( )
```

7.9.1.17 setup()

```
void setup ( )
```

7.9.1.18 stepperB()

```
AccelStepper stepperB (
    forwardstep1 ,
    backwardstep1 )
```

7.9.1.19 stepperT()

```
AccelStepper stepperT (
    forwardstep2 ,
    backwardstep2 )
```

7.9.2 Variable Documentation

7.9.2.1 acceleration

```
const float acceleration = 2000.0
```

Acceleration for the Stepper Motors.

7.9.2.2 AFMS

```
Adafruit_MotorShield AFMS = Adafruit_MotorShield()
```

7.9.2.3 collectionPattern

```
int collectionPattern = 1
```

[Collection](#) Pattern variable to choose between snake pattern XY Stage movement (default = 1) or Top down XY Stage movement (!= 1).

/** Arduino Setup. This only runs once.

7.9.2.4 columns

```
const int columns = 12
```

Amount of columns on well plate.

7.9.2.5 delayRPM

```
const int delayRPM = 255
```

Imports needed.

The Adafruit_MotorShieldV2 library and AccelStepper Library is used. Global Variables. Constains mostly motor values and digital pin values. RPM for the stepper motors.

7.9.2.6 limitSwitchB

```
const int limitSwitchB = 43
```

Digital Pin Number for Bottom Stepper Motor Limit Switch.

7.9.2.7 limitSwitchT

```
const int limitSwitchT = 41
```

Digital Pin Number for Top Stepper Motor Limit Switch.

7.9.2.8 locations

```
long locations[numLocations][2]
```

Three dimensional array for schematic of positions if MultiStepper is used.

7.9.2.9 maximumSpeed

```
const float maximumSpeed = 255.0
```

Max Speed of the Stepper Motors.

7.9.2.10 myMotorB

```
Adafruit_StepperMotor* myMotorB = AFMS.getStepper(stepsPerRevolution, 2)
```

7.9.2.11 myMotorT

```
Adafruit_StepperMotor* myMotorT = AFMS.getStepper(stepsPerRevolution, 1)
```

7.9.2.12 myPosition

```
int myPosition = 0
```

Position variable to keep track of which triplet well we are on.

7.9.2.13 numLocations

```
const int numLocations = ((rows * columns) / 3) + 1
```

Amount of triplet wells on well plate.

(Plus 1 for the eject location).

7.9.2.14 positions

```
long positions[2]
```

Array for movement needed for using both motors in MultiStepper.

7.9.2.15 rows

```
const int rows = 8
```

Amount of rows on well plate.

7.9.2.16 steppers

```
MultiStepper steppers
```

7.9.2.17 stepsPerRevolution

```
const int stepsPerRevolution = 200
```

Steps per revolution for the Stepper Motor.

7.9.2.18 stepsPerWell

```
const int stepsPerWell = 225
```

Amount of steps needed from center of well to the next.

7.9.2.19 stepsToEject

```
const int stepsToEject = 1000
```

Amount of steps to eject entire well plate.

7.9.2.20 TwoSwitchA

```
const int TwoSwitchA = 35
```

Digital Pin Number for TwoSwitchA.

7.9.2.21 TwoSwitchB

```
const int TwoSwitchB = 33
```

Digital Pin Number for TwoSwitchB.

7.9.2.22 TwoSwitchC

```
const int TwoSwitchC = 31
```

Digital Pin Number for TwoSwitchC.

Index

- `__init__`
 - `Collection.Collectador`, 26
 - `fakeSerial.Serial`, 45
 - `gui.App`, 14
 - `guiAutomaticPage`, 20
 - `guiManualPage`, 32
 - `guiSettingsPage`, 50
 - `guiWelcomePage`, 63
 - `Mswitch.MSwitch`, 41
 - `Pump.ThreePump`, 54
 - `TwoSwitch.TwoSwitch`, 60
 - `__str__`
 - `fakeSerial.Serial`, 45
- `abort_calibration`
 - `Pump.ThreePump`, 54
- `acceleration`
 - `ArduinoSnakePattern.cpp`, 72
 - `XYStage_LimitSwitch_TwoSwitch_Control.cpp`, 84
- `add_step`
 - `guiAutomaticPage`, 20
- `AFMS`
 - `ArduinoSnakePattern.cpp`, 72
 - `XYStage_LimitSwitch_TwoSwitch_Control.cpp`, 85
- `app`
 - `gui`, 11
- `appHighlightFont`
 - `gui.App`, 16
- `ArduinoSnakePattern.cpp`, 67
 - `acceleration`, 72
 - `AFMS`, 72
 - `backwardstep1`, 69
 - `backwardstep2`, 69
 - `columns`, 73
 - `delayRPM`, 73
 - `eject`, 69
 - `forwardstep1`, 69
 - `forwardstep2`, 69
 - `limitSwitchB`, 73
 - `limitSwitchT`, 73
 - `locations`, 73
 - `loop`, 70
 - `maximumSpeed`, 73
 - `moveDown`, 70
 - `moveGen`, 70
 - `moveLast`, 70
 - `moveLeft`, 71
 - `moveNext`, 71
 - `moveRight`, 71
 - `moveUp`, 71
 - `myMotorB`, 74
 - `myMotorT`, 74
 - `myPosition`, 74
 - `numLocations`, 74
 - `populateLocations`, 71
 - `positions`, 74
 - `reset`, 71
 - `rows`, 74
 - `set`, 72
 - `setup`, 72
 - `stepperB`, 72
 - `steppers`, 75
 - `stepperT`, 72
 - `stepsPerRevolution`, 75
 - `stepsPerWell`, 75
 - `stepsToEject`, 75
 - `TwoSwitchA`, 75
 - `TwoSwitchB`, 75
 - `TwoSwitchC`, 76
- `autocollectlabel`
 - `guiManualPage`, 36
- `autopagebutton`
 - `guiManualPage`, 36
 - `guiSettingsPage`, 50
- `backwardstep1`
 - `ArduinoSnakePattern.cpp`, 69
 - `XYStage_LimitSwitch_TwoSwitch_Control.cpp`, 81
- `backwardstep2`
 - `ArduinoSnakePattern.cpp`, 69
 - `XYStage_LimitSwitch_TwoSwitch_Control.cpp`, 81
- `baud`
 - `Collection.Collectador`, 29
 - `Mswitch.MSwitch`, 43
 - `Pump.ThreePump`, 58
 - `TwoSwitch.TwoSwitch`, 62
- `baudrate`
 - `fakeSerial.Serial`, 48
- `buttonFont`
 - `gui.App`, 16
- `bytesize`
 - `fakeSerial.Serial`, 48
- `calibrate`
 - `guiSettingsPage`, 50
 - `Pump.ThreePump`, 54
- `calibratebutton`
 - `guiSettingsPage`, 50
- `calibratecombobox`
 - `guiSettingsPage`, 51

- calibrationlabel
 - gui.SettingsPage, 51
- calibrationunit
 - Pump.ThreePump, 58
- calibrationvolume
 - Pump.ThreePump, 58
- channelcaliblabel
 - gui.SettingsPage, 51
- channels
 - gui.AutomaticPage, 22
 - gui.ManualPage, 36
 - gui.SettingsPage, 51
- channeltextlist
 - gui.ManualPage, 36
- channelvalues
 - gui.SettingsPage, 51
- check_all_manual
 - gui.ManualPage, 32
- check_all_manual_and_sampling
 - gui.ManualPage, 32
- check_channel
 - gui.AutomaticPage, 21
 - gui.ManualPage, 32
- check_input
 - gui, 10
- chop_return
 - Collection.Collectador, 26
 - Mswitch.MSwitch, 41
 - Pump.ThreePump, 55
 - TwoSwitch.TwoSwitch, 60
- clear_recipe
 - gui.AutomaticPage, 21
- clearbutton
 - gui.AutomaticPage, 22
- close
 - fakeSerial.Serial, 45
- coll_label
 - gui.WelcomePage, 64
- Collection, 9
- Collection.Collectador, 25
 - __init__, 26
 - baud, 29
 - chop_return, 26
 - connected, 29
 - currentPattern, 29
 - eject, 26
 - get_info, 26
 - last_site, 27
 - moveOneDown, 27
 - moveOneLeft, 27
 - moveOneRight, 27
 - moveOneUp, 27
 - next_site, 28
 - port, 30
 - position, 30
 - reset, 28
 - send, 28
 - ser, 30
 - serial_connect, 28
 - serialConnected, 30
 - setOrigin, 29
 - toggle_pattern, 29
 - uniqueID, 30
 - usingSnakePattern, 30
- Collection.py, 76
- CollectionControlLabel
 - gui.ManualPage, 36
- collectionPattern
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 85
- collectorframe
 - gui.WelcomePage, 64
- columns
 - ArduinoSnakePattern.cpp, 73
 - locationPopulate.cpp, 78
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 85
- combine_funcs
 - gui, 10
- connected
 - Collection.Collectador, 29
 - Mswitch.MSwitch, 43
 - Pump.ThreePump, 59
 - TwoSwitch.TwoSwitch, 62
- create_channels
 - gui.ManualPage, 33
- create_step
 - gui.AutomaticPage, 21
- currentPattern
 - Collection.Collectador, 29
- default_diameter_index
 - gui.App, 16
- delayRPM
 - ArduinoSnakePattern.cpp, 73
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 85
- delete_step
 - gui.AutomaticPage, 21
- deletestep
 - gui.AutomaticPage, 22
- deletestepcombo
 - gui.AutomaticPage, 22
- diameterlabel
 - gui.SettingsPage, 51
- diameterlist
 - gui.SettingsPage, 51
- diametervalues
 - gui.SettingsPage, 51
- directionheading
 - gui.SettingsPage, 52
- directionlabellist
 - gui.SettingsPage, 52
- directions
 - gui.SettingsPage, 52
- directionslist
 - gui.SettingsPage, 52
- done
 - gui.ManualPage, 33
- done_and_load

- gui.App, 14
- donebutton
 - gui.WelcomePage, 64
- eject
 - ArduinoSnakePattern.cpp, 69
 - Collection.Collectador, 26
 - gui.ManualPage, 33
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 81
- ejectbutton
 - guiAutomaticPage, 22
 - gui.ManualPage, 37
- entercalibratedbutton
 - gui.SettingsPage, 52
- fakeSerial, 9
- fakeSerial.py, 76
- fakeSerial.Serial, 44
 - __init__, 45
 - __str__, 45
 - baudrate, 48
 - bytesize, 48
 - close, 45
 - flush, 46
 - flushInput, 46
 - flushOutput, 46
 - isOpen, 46
 - name, 48
 - open, 46
 - parity, 48
 - port, 48
 - read, 47
 - readline, 47
 - rtscts, 48
 - stopbits, 48
 - timeout, 48
 - write, 47
 - xonxoff, 49
- file_load
 - guiAutomaticPage, 21
- file_save
 - guiAutomaticPage, 21
- flip
 - gui.WelcomePage, 64
- flowrateentrylist
 - gui.ManualPage, 37
- flush
 - fakeSerial.Serial, 46
- flushInput
 - fakeSerial.Serial, 46
- flushOutput
 - fakeSerial.Serial, 46
- FormatVolume
 - Pump.ThreePump, 55
- forwardstep1
 - ArduinoSnakePattern.cpp, 69
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 82
- forwardstep2
 - ArduinoSnakePattern.cpp, 69
- XYStage_LimitSwitch_TwoSwitch_Control.cpp, 82
- frames
 - gui.App, 17
- get_comport
 - gui.App, 14
- get_info
 - Collection.Collectador, 26
 - Mswitch.MSwitch, 42
- get_ports
 - gui.App, 15
 - Pump, 12
- getIdentifier
 - TwoSwitch.TwoSwitch, 60
- greaterlabelframe
 - gui.ManualPage, 37
- gui, 9
 - app, 11
 - check_input, 10
 - combine_funcs, 10
 - LARGE_FONT, 11
 - message_prompt, 10
- gui.App, 13
 - __init__, 14
 - appHighlightFont, 16
 - buttonFont, 16
 - default_diameter_index, 16
 - done_and_load, 14
 - frames, 17
 - get_comport, 14
 - get_ports, 15
 - hasColl, 17
 - hasMani, 17
 - hasPump, 17
 - max_flowrate_1, 17
 - max_flowrate_2, 17
 - max_flowrate_3, 17
 - max_flowrate_list, 18
 - message_window, 15
 - my2Switch, 18
 - myColl, 18
 - myMani, 18
 - myPump, 18
 - negative_run_status, 15
 - pageFont, 18
 - port2Switch, 18
 - portColl, 18
 - portMani, 19
 - portPump, 19
 - ports, 19
 - positive_run_status, 15
 - set_defaults, 16
 - should_be_running, 19
 - show_frame, 16
 - well_working_volume, 19
- guiAutomaticPage, 19
 - __init__, 20
 - add_step, 20
 - channels, 22

- check_channel, 21
- clear_recipe, 21
- clearbutton, 22
- create_step, 21
- delete_step, 21
- deletestep, 22
- deletestepcombo, 22
- ejectbutton, 22
- file_load, 21
- file_save, 21
- headinglabel, 23
- loadbutton, 23
- mainlabel, 23
- manualbutton, 23
- newstep, 23
- number_deleted, 23
- numberlist, 23
- reservoirs, 23
- run_recipe, 22
- runbutton, 24
- savebutton, 24
- settingsbutton, 24
- steplist, 24
- stepnumbers, 24
- todelete, 24
- twoSwitchValues, 24
- gui.ManualPage, 31
 - __init__, 32
 - autocollectlabel, 36
 - autopagebutton, 36
 - channels, 36
 - channeltextlist, 36
 - check_all_manual, 32
 - check_all_manual_and_sampling, 32
 - check_channel, 32
 - CollectionControlLabel, 36
 - create_channels, 33
 - done, 33
 - eject, 33
 - ejectbutton, 37
 - flowrateentrylist, 37
 - greaterlabelframe, 37
 - labelframelist, 37
 - lastwellbutton, 37
 - mainlabel, 37
 - moveOneDown, 33
 - moveOneLeft, 33
 - moveOneRight, 34
 - moveOneUp, 34
 - next, 34
 - nextwellbutton, 37
 - patternLabel, 37
 - positionLabel, 38
 - prev, 34
 - reservoircombobox, 38
 - reservoirlabel, 38
 - reservoirs, 38
 - reservoirvalues, 38
 - reset, 34
 - resetbutton, 38
 - sampleslabel, 38
 - samplingrate, 38
 - samplingrateentry, 39
 - set_and_start, 34
 - set_and_start_all, 35
 - setOrigin, 35
 - setOriginButton, 39
 - settingsbutton, 39
 - sites, 39
 - sitesentry, 39
 - start_and_collect, 35
 - startallbutton, 39
 - startandcollectbutton, 39
 - startbuttonlist, 39
 - stop, 35
 - stop_all, 35
 - stopallbutton, 40
 - stopbuttonlist, 40
 - toggle, 36
 - togglebutton, 40
 - toplables, 40
 - twoSwitchComboBoxList, 40
 - twoSwitchValues, 40
- gui.py, 76
- gui.SettingsPage, 49
 - __init__, 50
 - autopagebutton, 50
 - calibrate, 50
 - calibratebutton, 50
 - calibratecombobox, 51
 - calibrationlabel, 51
 - channelcaliblabel, 51
 - channels, 51
 - channelvalues, 51
 - diameterlabel, 51
 - diameterlist, 51
 - diametervvalues, 51
 - directionheading, 52
 - directionlabellist, 52
 - directions, 52
 - directionslist, 52
 - entercalibratedbutton, 52
 - mainlabel, 52
 - manualbutton, 52
 - measuredVolume, 52
 - measuredVolumeEntry, 53
 - reservoirs, 53
 - update, 50
 - updatebutton, 53
- gui.WelcomePage, 63
 - __init__, 63
 - coll_label, 64
 - collectorframe, 64
 - donebutton, 64
 - flip, 64
 - mainlabel, 64

- manilabel, 64
- perfusorframe, 64
- pumplabel, 65
- sublabelframe, 65
- hasColl
 - gui.App, 17
- hasMani
 - gui.App, 17
- hasPump
 - gui.App, 17
- headinglabel
 - guiAutomaticPage, 23
- isOn
 - Pump.ThreePump, 59
- isOpen
 - fakeSerial.Serial, 46
- labelframelist
 - gui.ManualPage, 37
- LARGE_FONT
 - gui, 11
- last_site
 - Collection.Collectador, 27
- lastwellbutton
 - gui.ManualPage, 37
- limitSwitchB
 - ArduinoSnakePattern.cpp, 73
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 85
- limitSwitchT
 - ArduinoSnakePattern.cpp, 73
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 85
- loadbutton
 - guiAutomaticPage, 23
- locationPopulate.cpp, 77
 - columns, 78
 - main, 77
 - numLocations, 78
 - populateLocations, 77
 - rows, 78
 - stepsPerWell, 78
 - stepsToEject, 78
- locations
 - ArduinoSnakePattern.cpp, 73
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 86
- loop
 - ArduinoSnakePattern.cpp, 70
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 82
- main
 - locationPopulate.cpp, 77
- mainlabel
 - guiAutomaticPage, 23
 - gui.ManualPage, 37
 - gui.SettingsPage, 52
 - gui.WelcomePage, 64
- manilabel
 - gui.WelcomePage, 64
- manualbutton
 - guiAutomaticPage, 23
 - gui.SettingsPage, 52
- max_flowrate_1
 - gui.App, 17
- max_flowrate_2
 - gui.App, 17
- max_flowrate_3
 - gui.App, 17
- max_flowrate_list
 - gui.App, 18
- maximumSpeed
 - ArduinoSnakePattern.cpp, 73
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 86
- measuredVolume
 - gui.SettingsPage, 52
- measuredVolumeEntry
 - gui.SettingsPage, 53
- message_prompt
 - gui, 10
- message_window
 - gui.App, 15
- moveDown
 - ArduinoSnakePattern.cpp, 70
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 82
- moveGen
 - ArduinoSnakePattern.cpp, 70
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 82
- moveLast
 - ArduinoSnakePattern.cpp, 70
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 83
- moveLeft
 - ArduinoSnakePattern.cpp, 71
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 83
- moveNext
 - ArduinoSnakePattern.cpp, 71
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 83
- moveOneDown
 - Collection.Collectador, 27
 - gui.ManualPage, 33
- moveOneLeft
 - Collection.Collectador, 27
 - gui.ManualPage, 33
- moveOneRight
 - Collection.Collectador, 27
 - gui.ManualPage, 34
- moveOneUp
 - Collection.Collectador, 27
 - gui.ManualPage, 34
- moveRight
 - ArduinoSnakePattern.cpp, 71
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 83
- moveUp
 - ArduinoSnakePattern.cpp, 71
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 83
- Mswitch, 11
- Mswitch.MSwitch, 41
 - __init__, 41

- baud, 43
- chop_return, 41
- connected, 43
- get_info, 42
- port, 43
- res, 43
- send, 42
- ser, 43
- serial_connect, 42
- serialConnected, 43
- set_reservoir, 42
- uniqueID, 44
- Mswitch.py, 79
- my2Switch
 - gui.App, 18
- myColl
 - gui.App, 18
- myMani
 - gui.App, 18
- myMotorB
 - ArduinoSnakePattern.cpp, 74
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 86
- myMotorT
 - ArduinoSnakePattern.cpp, 74
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 86
- myPosition
 - ArduinoSnakePattern.cpp, 74
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 86
- myPump
 - gui.App, 18
- name
 - fakeSerial.Serial, 48
- negative_run_status
 - gui.App, 15
- newstep
 - gui.AutomaticPage, 23
- next
 - gui.ManualPage, 34
- next_site
 - Collection.Collectador, 28
- nextwellbutton
 - gui.ManualPage, 37
- number_deleted
 - gui.AutomaticPage, 23
- numberlist
 - gui.AutomaticPage, 23
- numLocations
 - ArduinoSnakePattern.cpp, 74
 - locationPopulate.cpp, 78
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 86
- open
 - fakeSerial.Serial, 46
- pageFont
 - gui.App, 18
- parity
 - fakeSerial.Serial, 48
- patternLabel
 - gui.ManualPage, 37
- perfusorframe
 - gui.WelcomePage, 64
- populateLocations
 - ArduinoSnakePattern.cpp, 71
 - locationPopulate.cpp, 77
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 83
- port
 - Collection.Collectador, 30
 - fakeSerial.Serial, 48
 - Mswitch.MSwitch, 43
 - Pump.ThreePump, 59
 - TwoSwitch.TwoSwitch, 62
- port2Switch
 - gui.App, 18
- portColl
 - gui.App, 18
- portMani
 - gui.App, 19
- portPump
 - gui.App, 19
- ports
 - gui.App, 19
- position
 - Collection.Collectador, 30
- positionLabel
 - gui.ManualPage, 38
- positions
 - ArduinoSnakePattern.cpp, 74
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 87
- positive_run_status
 - gui.App, 15
- prev
 - gui.ManualPage, 34
- Pump, 11
 - get_ports, 12
- Pump.py, 79
- Pump.ThreePump, 53
 - __init__, 54
 - abort_calibration, 54
 - baud, 58
 - calibrate, 54
 - calibrationunit, 58
 - calibrationvolume, 58
 - chop_return, 55
 - connected, 59
 - FormatVolume, 55
 - isOn, 59
 - port, 59
 - send, 55
 - send_return, 55
 - ser, 59
 - serial_connect, 56
 - serialConnected, 59
 - setDefaults, 56
 - setDir, 56
 - setFlow, 56

- setMeasuredVolume, 57
 - setTargetCalibrationVolume, 57
 - start, 57
 - start_all, 57
 - stop, 58
 - stop_all, 58
- pumplabel
 - gui.WelcomePage, 65
- read
 - fakeSerial.Serial, 47
- readline
 - fakeSerial.Serial, 47
- res
 - Mswitch.MSwitch, 43
- reservoircombobox
 - gui.ManualPage, 38
- reservoirlabel
 - gui.ManualPage, 38
- reservoirs
 - guiAutomaticPage, 23
 - gui.ManualPage, 38
 - gui.SettingsPage, 53
- reservoirvalues
 - gui.ManualPage, 38
- reset
 - ArduinoSnakePattern.cpp, 71
 - Collection.Collectador, 28
 - gui.ManualPage, 34
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 84
- resetbutton
 - gui.ManualPage, 38
- rows
 - ArduinoSnakePattern.cpp, 74
 - locationPopulate.cpp, 78
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 87
- rtscts
 - fakeSerial.Serial, 48
- run_recipe
 - guiAutomaticPage, 22
- runbutton
 - guiAutomaticPage, 24
- sampleslabel
 - gui.ManualPage, 38
- samplingrate
 - gui.ManualPage, 38
- samplingrateentry
 - gui.ManualPage, 39
- savebutton
 - guiAutomaticPage, 24
- send
 - Collection.Collectador, 28
 - Mswitch.MSwitch, 42
 - Pump.ThreePump, 55
 - TwoSwitch.TwoSwitch, 61
- send_return
 - Pump.ThreePump, 55
- ser
 - Collection.Collectador, 30
 - Mswitch.MSwitch, 43
 - Pump.ThreePump, 59
 - TwoSwitch.TwoSwitch, 62
- serial_connect
 - Collection.Collectador, 28
 - Mswitch.MSwitch, 42
 - Pump.ThreePump, 56
 - TwoSwitch.TwoSwitch, 61
- serialConnected
 - Collection.Collectador, 30
 - Mswitch.MSwitch, 43
 - Pump.ThreePump, 59
 - TwoSwitch.TwoSwitch, 62
- set
 - ArduinoSnakePattern.cpp, 72
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 84
- set_and_start
 - gui.ManualPage, 34
- set_and_start_all
 - gui.ManualPage, 35
- set_defaults
 - gui.App, 16
- set_reservoir
 - Mswitch.MSwitch, 42
- setCollect
 - TwoSwitch.TwoSwitch, 61
- setDefault
 - Pump.ThreePump, 56
- setDir
 - Pump.ThreePump, 56
- setFlow
 - Pump.ThreePump, 56
- setMeasuredVolume
 - Pump.ThreePump, 57
- setOrigin
 - Collection.Collectador, 29
 - gui.ManualPage, 35
- setOriginButton
 - gui.ManualPage, 39
- setRecirculate
 - TwoSwitch.TwoSwitch, 61
- setTargetCalibrationVolume
 - Pump.ThreePump, 57
- settingsbutton
 - guiAutomaticPage, 24
 - gui.ManualPage, 39
- setup
 - ArduinoSnakePattern.cpp, 72
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 84
- should_be_running
 - gui.App, 19
- show_frame
 - gui.App, 16
- sites
 - gui.ManualPage, 39
- sitesentry
 - gui.ManualPage, 39

- start
 - Pump.ThreePump, 57
- start_all
 - Pump.ThreePump, 57
- start_and_collect
 - gui.ManualPage, 35
- startallbutton
 - gui.ManualPage, 39
- startandcollectbutton
 - gui.ManualPage, 39
- startbuttonlist
 - gui.ManualPage, 39
- steplist
 - gui.AutomaticPage, 24
- stepnumbers
 - gui.AutomaticPage, 24
- stepperB
 - ArduinoSnakePattern.cpp, 72
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 84
- steppers
 - ArduinoSnakePattern.cpp, 75
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 87
- stepperT
 - ArduinoSnakePattern.cpp, 72
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 84
- stepsPerRevolution
 - ArduinoSnakePattern.cpp, 75
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 87
- stepsPerWell
 - ArduinoSnakePattern.cpp, 75
 - locationPopulate.cpp, 78
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 87
- stepsToEject
 - ArduinoSnakePattern.cpp, 75
 - locationPopulate.cpp, 78
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 87
- stop
 - gui.ManualPage, 35
 - Pump.ThreePump, 58
- stop_all
 - gui.ManualPage, 35
 - Pump.ThreePump, 58
- stopallbutton
 - gui.ManualPage, 40
- stopbits
 - fakeSerial.Serial, 48
- stopbuttonlist
 - gui.ManualPage, 40
- sublabelframe
 - gui.WelcomePage, 65
- timeout
 - fakeSerial.Serial, 48
- todelete
 - gui.AutomaticPage, 24
- toggle
 - gui.ManualPage, 36
- toggle_pattern
 - Collection.Collectador, 29
- togglebutton
 - gui.ManualPage, 40
- toplabels
 - gui.ManualPage, 40
- TwoSwitch, 12
- TwoSwitch.py, 79
- TwoSwitch.TwoSwitch, 59
 - __init__, 60
 - baud, 62
 - chop_return, 60
 - connected, 62
 - getIdentifier, 60
 - port, 62
 - send, 61
 - ser, 62
 - serial_connect, 61
 - serialConnected, 62
 - setCollect, 61
 - setRecirculate, 61
 - uniqueID, 62
 - willRecirculate, 63
- TwoSwitchA
 - ArduinoSnakePattern.cpp, 75
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 88
- TwoSwitchB
 - ArduinoSnakePattern.cpp, 75
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 88
- TwoSwitchC
 - ArduinoSnakePattern.cpp, 76
 - XYStage_LimitSwitch_TwoSwitch_Control.cpp, 88
- twoSwitchComboBoxList
 - gui.ManualPage, 40
- twoSwitchValues
 - gui.AutomaticPage, 24
 - gui.ManualPage, 40
- uniqueID
 - Collection.Collectador, 30
 - Mswitch.MSwitch, 44
 - TwoSwitch.TwoSwitch, 62
- update
 - gui.SettingsPage, 50
- updatebutton
 - gui.SettingsPage, 53
- usingSnakePattern
 - Collection.Collectador, 30
- well_working_volume
 - gui.App, 19
- willRecirculate
 - TwoSwitch.TwoSwitch, 63
- write
 - fakeSerial.Serial, 47
- xonxoff
 - fakeSerial.Serial, 49
- XYStage_LimitSwitch_TwoSwitch_Control.cpp, 79
 - acceleration, 84
 - AFMS, 85

- backwardstep1, [81](#)
- backwardstep2, [81](#)
- collectionPattern, [85](#)
- columns, [85](#)
- delayRPM, [85](#)
- eject, [81](#)
- forwardstep1, [82](#)
- forwardstep2, [82](#)
- limitSwitchB, [85](#)
- limitSwitchT, [85](#)
- locations, [86](#)
- loop, [82](#)
- maximumSpeed, [86](#)
- moveDown, [82](#)
- moveGen, [82](#)
- moveLast, [83](#)
- moveLeft, [83](#)
- moveNext, [83](#)
- moveRight, [83](#)
- moveUp, [83](#)
- myMotorB, [86](#)
- myMotorT, [86](#)
- myPosition, [86](#)
- numLocations, [86](#)
- populateLocations, [83](#)
- positions, [87](#)
- reset, [84](#)
- rows, [87](#)
- set, [84](#)
- setup, [84](#)
- stepperB, [84](#)
- steppers, [87](#)
- stepperT, [84](#)
- stepsPerRevolution, [87](#)
- stepsPerWell, [87](#)
- stepsToEject, [87](#)
- TwoSwitchA, [88](#)
- TwoSwitchB, [88](#)
- TwoSwitchC, [88](#)